

What doctors are saying:

"In a single day I was able to do 11 mock-up veneers that were not part of the original schedule, resulting in several new cases being accepted. I wouldn't want to work without them now."

-Dr. Chad Wagener, DDS

"Terrific tool to quickly and easily create beautiful anterior restorations."

-Dr. Gary M. Radz, DDS

"The simplicity of Uveneer is absolutely remarkable. Why didn't I think of this?"

-Dr. John C. Comisi, DDS, MAGD


"Uveneer makes the practice of the anterior esthetic dentistry easier, faster, and better."

-Dr. George Freedman, DDS



UVKV3 Uveneer Kit
16 x Medium upper and lower arch templates
16 x Large upper and lower arch templates

Medium and Large templates provide 2 central incisors, 2 lateral incisors, 2 canines, and 2 premolar templates for both the upper and lower arches.



UVKEV1 Uveneer Extra Kit
6 x Extra Large upper anterior templates
6 x Large upper anterior templates
6 x Medium upper anterior templates
6 x Square upper anterior templates

UVKEXLSQV1 Uveneer Extra XL & SQ Kit
6 x Extra Large upper anterior templates
6 x Square upper anterior templates

UVKELMV1 Uveneer Extra L & M Kit
6 x Large upper anterior templates
6 x Medium upper anterior templates



Follow Us!

Join the Uveneer template Facebook group to share cases, learn techniques, and ask Dr. Sigal Jacobson questions!

[facebook.com/groups/uveneer](https://www.facebook.com/groups/uveneer)

Follow the Uveneer template brand for news, updates, and videos.

[facebook.com/uveneerdentist](https://www.facebook.com/uveneerdentist)



SCAN TO VIEW PRICING



DIRECT COMPOSITE TEMPLATE SYSTEMS



Simple and Versatile

The Complete Uveneer™ Template Family

The Uveneer and Uveneer Extra™ template systems are unique, minimally invasive, and create beautiful direct composite veneers with predictable shape and symmetry. They mimic natural dentition and are designed to create high-quality anterior restorations in one visit. They can also be used for cosmetic mock-ups and shade selection, as well as temporaries during porcelain veneer creation.



Uveneer Templates Promise a Beautiful Symmetry



A Way to Provide Unique and Natural Smiles



The original Uveneer template kit contains upper and lower anterior templates in Medium and Large which:

- Create predictable, reproducible composite restorations
- Prevent the oxygen inhibition layer during curing, resulting in a hard and glossy surface
- Allow light to pass through the template to the composite for effective curing
- Work with any preferred composite
- Release easily from cured composite resin
- Facilitate application on individual or multiple teeth
- Are autoclavable* and reusable, making them a cost-effective choice

*Black base and transparent lids are NOT autoclavable. Only templates are autoclavable. See instructions for use for more information.

Uveneer Extra templates provide the option of four new upper anterior designs: Extra Large, Large, Medium, and Square. To complement the original Uveneer templates, they also:

- Create natural-looking contours on composite restorations
- Require minimal adjusting or polishing, saving time
- Mimic real tooth anatomy for authentic, natural results
- Are molded based on scanned teeth from Dr. Jan Hajto's work
- Add to the unique potential of the patient's smile
- Possess all of the trusted qualities of the original Uveneer templates, but with a new design



L - Large Upper

L - Large Lower

M - Medium Upper

M - Medium Lower

XL - Extra Large

L - Large

M - Medium

SQ - Square