

## What doctors are saying:

"In a single day I was able to do 11 mock-up veneers that were not part of the original schedule, resulting in several new cases being accepted. I wouldn't want to work without them now."

**-Dr. Chad Wagener, DDS**

"Terrific tool to quickly and easily create beautiful anterior restorations."

**-Dr. Gary M. Radz, DDS**

"The simplicity of Uveneer is absolutely remarkable. Why didn't I think of this?"

**-Dr. John C. Comisi, DDS, MAGD**


"Uveneer makes the practice of the anterior esthetic dentistry easier, faster, and better."

**-Dr. George Freedman, DDS**



**UVKV3 Uveneer Kit**  
16 x Medium upper and lower arch templates  
16 x Large upper and lower arch templates

Medium and Large templates provide 2 central incisors, 2 lateral incisors, 2 canines, and 2 premolar templates for both the upper and lower arches.



**UVKEV1 Uveneer Extra Kit**  
6 x Extra Large upper anterior templates  
6 x Large upper anterior templates  
6 x Medium upper anterior templates  
6 x Square upper anterior templates

**UVKEXLSQV1 Uveneer Extra XL & SQ Kit**  
6 x Extra Large upper anterior templates  
6 x Square upper anterior templates

**UVKELMV1 Uveneer Extra L & M Kit**  
6 x Large upper anterior templates  
6 x Medium upper anterior templates

**Uveneer™**  
dental art made simple

**Uveneer™**  
EXTRA

DIRECT COMPOSITE TEMPLATE SYSTEMS



Simple and Versatile

## Follow Us!

Join the Uveneer template Facebook group to share cases, learn techniques, and ask Dr. Sigal Jacobson questions!

[facebook.com/groups/uveneer](https://www.facebook.com/groups/uveneer)

Follow the Uveneer template brand for news, updates, and videos.

[facebook.com/uveneerdentist](https://www.facebook.com/uveneerdentist)



# The Complete Uveneer™ Template Family

The Uveneer and Uveneer Extra™ template systems are unique, minimally invasive, and create beautiful direct composite veneers with predictable shape and symmetry. They mimic natural dentition and are designed to create high-quality anterior restorations in one visit. They can also be used for cosmetic mock-ups and shade selection, as well as temporaries during porcelain veneer creation.



## Uveneer Templates Promise a Beautiful Symmetry



## A Way to Provide Unique and Natural Smiles



**The original Uveneer template kit contains upper and lower anterior templates in Medium and Large which:**

- Create predictable, reproducible composite restorations
- Prevent the oxygen inhibition layer during curing, resulting in a hard and glossy surface
- Allow light to pass through the template to the composite for effective curing
- Work with any preferred composite
- Release easily from cured composite resin
- Facilitate application on individual or multiple teeth
- Are autoclavable\* and reusable, making them a cost-effective choice

\*Black base and transparent lids are NOT autoclavable. Only templates are autoclavable. See instructions for use for more information.

**Uveneer Extra templates provide the option of four new upper anterior designs: Extra Large, Large, Medium, and Square. To complement the original Uveneer templates, they also:**

- Create natural-looking contours on composite restorations
- Require minimal adjusting or polishing, saving time
- Mimic real tooth anatomy for authentic, natural results
- Are molded based on scanned teeth from Dr. Jan Hajto's work
- Add to the unique potential of the patient's smile
- Possess all of the trusted qualities of the original Uveneer templates, but with a new design



L - Large Upper

L - Large Lower

M - Medium Upper

M - Medium Lower

XL - Extra Large

L - Large

M - Medium

SQ - Square