



50%
Bigger
Lens

VALO[®]

GRAND

ECLIPSING THE COMPETITION.



ULTRADENT
PRODUCTS, INC.

800.552.5512 ULTRADENT.COM



ALL OF THE BENEFITS OF THE REVOLUTIONARY,
AWARD-WINNING VALO® WITH A 50% BIGGER LENS.



- Curing button on the underside for intuitive operation
- Optimally collimated beam delivers consistent, uniform power
- Three curing modes—Standard, High Power Plus, and Xtra Power—accommodate your preferences
- Unibody construction yields a durable, slim, ergonomic shape and allows unprecedented access to all restoration sites
- Highly efficient LEDs keep the wand body cool to the touch
- Battery-operated, cordless wand design provides optimal convenience and flexibility
- Rechargeable batteries are safe and inexpensive



**LARGER
HEAD SIZE**



**ADDITIONAL
CURING
BUTTON**



**BROADBAND
LED**



**DURABLE
CONSTRUCTION**



**ACCESS
ANGLE**

5972 - VALO Grand Kit

- 1 x VALO Grand LED curing light
- 4 x Batteries (rechargeable)
- 1 x Battery charger
- 1 x Power supply with universal plugs for battery charger
- 1 x VALO Grand light shield
- 1 x Handpiece bracket holder
- 1 x Sample pack of barrier sleeves

ACCESSORIES:

- 4666** - VALO Grand Barrier Sleeves 100pk
- 5961** - VALO Cordless Charging Unit Power Supply
- 5962** - VALO Cordless Battery Charging Unit
- 5963** - VALO Cordless Rechargeable Batteries 2pk
- 1667** - Surface Mounting Bracket
- 3604** - VALO Grand Cordless Light Shield* 1pk

*Reusable. Should be disinfected using an intermediate-level disinfectant.

800.552.5512 ULTRADENT.COM



© 2019 Ultradent Products, Inc. All rights reserved.
1004596.4 101118

ULTRADENT
PRODUCTS, INC.

Improving Oral Health Globally