

BFC3TM
POWERED IMPRESSION GUN



Meet the

Motorized Portable Impression Gun

Better Faster Comfort, Control, Compatibility—The BFC3

The Ultradent BFC3 impression gun makes easy work of impression taking. Our BFC3 gun is the first cordless, handheld, motorized impression gun that dispenses standard 25 ml, 50 ml, and 75 ml 1:1 or 2:1 automix cartridges.

BFC3™ POWERED IMPRESSION GUN

Try out the BFC3 gun and auto-dispense impression materials chairside. You won't be constrained by power cords or heavy, tabletop automated dispensers. Dispense and prepare impressions conveniently and efficiently.

Comfort

Dispense automix material with a press of a button! It can fill a full-arch tray with heavy-body material in less than 10 seconds. No longer will you need to repeatedly squeeze a large spring-loaded trigger to express impression material into a tray.

Control

An anti-drip feature momentarily reverses the plunger when the dispensing trigger is released. Added to that is an auto-stop feature with overload protection when a cartridge has been emptied.

The adjustable flow rate on the BFC3 gun optimizes the expression of various impression material bodies.

Compatibility

Use your existing products, and don't worry about buying any special cartridges. This gun is convenient and compatible with the system you have in place at your office.



Additional Features:

- Ergonomic design similar to automix guns
- Includes convenient countertop stand for easy storage and access
- 1-year warranty
- On a fully charged battery, the Ultradent BFC3 Powered Impression Gun can dispense about 60 full 50 ml heavy body VPS cartridges
- Low battery indicator light.
(Ultradent BFC3 Powered Impression Gun will dispense up to 10 full cartridges after the low battery LED illuminates.)
- Battery charger features full-charge indicator LED and overcharge protection circuitry
- Ultradent BFC3 battery is rated for 500 recharges (about 30,000 cartridges) before a replacement will be needed



5565 - BFC3 Impression Material Gun Kit

- 1 x Ultradent BFC3 Powered Impression Gun
- 1 x Countertop Stand
- 1 x 50 ml/25 ml Plunger
- 1 x 75 ml Plunger
- 1 x Battery Charger



5568 - BFC3 Replacement Plunger Set

- 1 x 50 ml/25 ml Plunger
- 1 x 75 ml Plunger



5569 - BFC3 Replacement Battery Charger

- 1 x Battery Charger

Product Specifications

BATTERY:	Rechargeable 11.1V lithium-ion
BATTERY CHARGER:	12.6V/500mA with LED indicator
CHARGING TIME:	2 1/2 hours (full charge)
FLOW RATE:	Adjustable from 0.6 to 1.7ml/sec.
DEVICE DIMENSIONS:	4"W x 4"D x 7"H (overall)
DEVICE WEIGHT:	28.1 oz. (without cartridge)

Call your Ultradent Sales Representative to order your own! Or shop at ultradent.com



SCAN
TO VIEW
PRICING



© 2023 Ultradent Products, Inc. All rights reserved.
1011779AR01 012523

800.552.5512 | ULTRADENT.COM

