

# Perfect Delivery. Optimal Performance. Easy Removal.

simplifying



endodontics

## UltraCal™ XS Paste and Ultradent™ Citric Acid 20% Solution

**UltraCal XS calcium hydroxide paste** (pH 12.5) is easily delivered with the NaviTip™ tip exactly where it is needed in the canal. Its aqueous consistency and radiopaque properties make it ideal as an endodontic treatment dressing for two-appointment RCTs.

When it comes time to remove UltraCal XS paste from the canal, look no further than **Ultradent Citric Acid Solution**, delivered with the NaviTip™ FX™ tip. Citric Acid solution dissolves calcium hydroxide, and the small fibers attached to the NaviTip FX tip easily scrub the walls of the canal, removing any remaining calcium hydroxide as well as the smear layer.



# UltraCal™ XS

## UltraCal™ XS

30%–35% CALCIUM HYDROXIDE PASTE



**1027 - UltraCal XS Syringe Kit**  
 4 x 1.2 ml syringes  
 5 x Each 29 ga NaviTip tips  
 17 mm, 21 mm, 25 mm, and 27 mm



**606 - UltraCal XS Syringe 4pk**  
**5117 - UltraCal XS Syringe 20pk**  
 1.2 ml syringes



**329 - Citric Acid IndiSpense Syringe 1pk**  
 30 ml syringe


## NaviTip™ Tips



LOK-TITE™	mm	20pk	50pk
29 ga NaviTip	27	5115	1377
29 ga NaviTip	25	5114	1376
29 ga NaviTip	21	5113	1374
29 ga NaviTip	17	5112	1378
29 ga–29 ga NaviTips	27–17	5116	1379

LOK-TITE™	mm	20pk	50pk
30 ga NaviTip	27	1354	1424
30 ga NaviTip	25	1250	1423
30 ga NaviTip	21	1349	1422
30 ga NaviTip	17	1249	1421
30 ga–30 ga NaviTips	27–17	1351	3319

## NaviTip™ FX™ Tip



LOK-TITE™	mm	20pk
30 ga NaviTip FX	17	1452
30 ga NaviTip FX	25	1454

