

ultrapro™ Tx
sweep™



REDUCE SPLATTER UP TO 95%¹




The Ultrapro Tx Sweep disposable prophy angle features a new, innovative brush guard that prevents saliva and paste from collecting and roping around the cup, which helps prevent splatter. Five out of six hygienists experienced a significant reduction in splatter when using the Sweep prophy angle.²

1. Data on file.
2. Data on file.


The Ultrapro Tx Sweep disposable prophylaxis angle has an innovative brush guard with flexible bristles that efficiently keep saliva from collecting and roping around the cup. The cup also sweeps the prophylaxis paste back toward the tooth, so the paste is not wasted and lost in the saliva.

- Brush helps prevent accumulation of saliva and paste on the outside of the cup
- Sweeps paste back toward the tooth, allowing the clinician to continue working and make fewer stops to refill the cup with more prophylaxis paste
- Advanced internal blade design
- Outer ridges for interproximal cleaning
- A 20% shorter head and 25% slimmer neck design gives better access and improved visibility
- Optimal flare
- Comfortable, ergonomic shape
- Not made with rubber latex
- Available in soft or firm cup design




Splatter Comparison Testing

Sweep prophylaxis angle



Traditional prophylaxis angle



sweep
(Available in soft or firm cup)

Ultrapro Tx Sweep Prophylaxis Angles

8357 Soft, Green 144pk

8359 Soft, Green 500pk



8358 Firm, Blue 144pk

8360 Firm, Blue 500pk

