



**8502 - Ultrapro Tx Air Kit**  
 1 x Motor  
 3 x Nose cones  
 1 x Universal lubricating oil



**8506 - Ultrapro Tx Air Motor 1pk**



**8503 - Ultrapro Tx Air Green Nose Cone 1pk**  
**8504 - Ultrapro Tx Air Pink Nose Cone 1pk**  
**8505 - Ultrapro Tx Air Blue Nose Cone 1pk**

**8316 - Universal Lubricating Oil - 0.5 oz Bottle 1pk**



**5533 - Ultrapro Tx Cordless Prophy Foot Control Replacement Part**  
 1 x Foot control



**5532 - Ultrapro Tx Cordless Prophy O-Ring Replacement Part**  
 6 x Black sheath O-rings  
 2 x Green motor O-rings



**5275 - Ultrapro Tx Cordless Prophy Motor Kit**  
 1 x Cordless prophy motor  
 3 x Sheaths  
 1 x Foot control  
 1 x Charger base  
 2 x Wall chargers  
 4 x Prophy angles  
 100 x Barrier sleeves  
 6 x Black sheath O-rings  
 2 x Green motor O-rings

**5276 - Ultrapro Tx Cordless Prophy Sheath Kit**  
 1 x Sheath  
 4 x Black sheath O-rings  
 2 x Green motor O-rings



**5531 - Ultrapro Tx Cordless Prophy Wall Charger Replacement Part**  
 1 x Wall Charger



**5499 - Ultrapro Tx Cordless Prophy Motor Barrier Sleeves 500pk**  
 500 x Barrier sleeves



**5534 - Ultrapro Tx Cordless Prophy Charger Base Replacement Part**  
 1 x Charger base



# TRUE INNOVATION FOR EVERY POLISHING NEED



It's a complete line for your practice. Our Ultrapro™ Tx Air handpiece and the new Ultrapro™ Tx Cordless handpiece are designed to be light, ergonomic, and comfortable. Plus, they work with nearly any DPA for high versatility. Our Ultrapro™ Tx Sweep™ DPA significantly reduces splatter—a truly innovative solution as part of the whole Ultrapro Tx product line.



SCAN TO VIEW PRICING

# ultrapro™ Tx air

Designed with comfort in mind, both you and your patient will enjoy the smooth polishing experience. Lightweight and ergonomic, the Ultrapro Tx Air handpiece will be an efficient and effortless tool you'll love to use.

- Lightweight—(2.2 oz) to reduce hand and arm fatigue<sup>1</sup>
- Ergonomically shaped and low vibration for maximum comfort
- 360° swivel for easy rotation
- Soft-start motor to minimize splatter
- Maximum speed of 4,000 RPM
- Universal E-type motor for multifunctionality
- Compatible with nearly all disposable prophylaxis angles



blue  
green  
pink



1. realityesthetics.com.

# ultrapro™ Tx cordless

Meet the new Ultrapro™ Tx Cordless handpiece—the hygiene handpiece that gives you the power of a corded handpiece with the unlimited movement of cordless. The cordless design means you have the ability to approach a tooth from any angle with no cord drag.

- Cordless design allows for unlimited movement
- Compatible with most DPAs for maximum usability
- A single wireless charge can last an entire day
- Automatic shutoff after five minutes
- Working output of up to 4,000 RPM
- Wireless Bluetooth foot control
- Lightweight and ergonomic
- 2-year warranty



# ultrapro™ Tx disposable prophylaxis angles

Ultrapro™ Tx disposable prophylaxis angles are innovative, effective, and comfortable to use. They feature smooth, quiet gears and an ergonomic design so that both you and your patient will have a comfortable experience. With an advanced cup design both inside and out, the Extra, Skini, and Sweep angles are built to clean better than ever. The Ultrapro Tx Sweep disposable prophylaxis angle features a new, innovative brush guard that prevents saliva and paste from collecting and roping around the cup to help reduce splatter up to 95%.<sup>2</sup>

- Innovative cup design for reduced splatter and efficient cleaning
- Outer ridges for improved interproximal cleaning (Extra, Skini, and Sweep angles)
- Designed for better access and improved visibility, with 20% shorter head and 25% slimmer neck (Skini and Sweep angles)
- Available in right-angle or contra-angle design (Extra and Skini angles)
- All designs feature optimal flare
- Ergonomic design for comfort
- Not made with natural rubber latex
- Available in soft or firm cup designs
- Skini angle is sterile to minimize risk of infection
- Sweep angle features an innovative brush design that significantly reduces splatter

### DPA Head Comparison



# ultrapro™ Tx prophy paste

Ultrapro™ Tx prophy paste and Ultrapro™ Tx Pure prophy paste perfectly round out all your cleaning needs. Low splatter and easy to rinse, they provide a comfortable cleaning experience. With a variety of flavors to choose from, Ultrapro Tx prophy pastes will impress your patients with their great taste and leave them looking forward to their next visit.




### Ultrapro Tx 2 g 200pk

| Flavor                      | Fine | Medium | Coarse |
|-----------------------------|------|--------|--------|
| Cool Mint                   | 8309 | 8310   | 8311   |
| Walterberry                 | 8320 | 8312   | 8313   |
| Bubble Gum dye-free         | 8321 | 8314   | 8315   |
| Orange Cream                | 8322 | 8323   | 8324   |
| Variety W-100, OD-50, BG-50 | 8326 | —      | —      |
| Variety W-100, OD-50, CM-50 | —    | 8327   | —      |
| PURE                        | —    | 8325   | —      |




1. Moravec, L. Habits for a healthy practice. *Dimen Dent Hyg.* 2012;10(11):42–44. 2. Data on file.


### 8356 - Ultrapro Tx Tapered Brush 100pk



| original | 144pk | 500pk |
|----------|-------|-------|
| Soft     | 8308  | 8318  |
| Firm     | 8307  | 8317  |



| extra | 144pk | 500pk |
|-------|-------|-------|
| Soft  | 8340  | 8342  |
| Firm  | 8341  | 8343  |



| extra CONTRA | 144pk | 500pk |
|--------------|-------|-------|
| Soft         | 8344  | 8346  |
| Firm         | 8345  | 8347  |



| sterile skini | 144pk | 500pk |
|---------------|-------|-------|
| Soft          | 8348  | 8350  |
| Firm          | 8349  | 8351  |



| sterile skini CONTRA | 144pk | 500pk |
|----------------------|-------|-------|
| Soft                 | 8352  | 8354  |
| Firm                 | 8353  | 8355  |



| sweep | 144pk | 500pk |
|-------|-------|-------|
| Soft  | 8357  | 8359  |
| Firm  | 8358  | 8360  |

