




8502 - Ultrapro Tx Air Kit  
1 x Motor  
3 x Nose cones  
1 x Universal lubricating oil



8506 - Ultrapro Tx Air Motor 1pk



8503 - Ultrapro Tx Air Green Nose Cone 1pk  
8504 - Ultrapro Tx Air Pink Nose Cone 1pk  
8505 - Ultrapro Tx Air Blue Nose Cone 1pk



8316 - Universal Lubricating Oil - 0.5 oz Bottle 1pk



COOL MINT   WALTERBERRY™   BUBBLE GUM   ORANGE DREAM-SICLE   PURE


Ultrapro Tx 2 g 200pk

Flavor	Fine	Medium	Coarse
Cool Mint	8309	8310	8311
Walterberry	8320	8312	8313
Bubble Gum dye-free	8321	8314	8315
Orange Dreamsicle	8322	8323	8324
Variety W-100, OD-50, BG-50	8326	—	—
Variety W-100, OD-50, CM-50	—	8327	—
PURE	—	8325	—

8356 - Ultrapro Tx Tapered Brush 100pk



original	144pk	500pk
Soft	8308	8318
Firm	8307	8317



extra	144pk	500pk
Soft	8340	8342
Firm	8341	8343



extra CONTRA	144pk	500pk
Soft	8344	8346
Firm	8345	8347



sterile R skini	144pk	500pk
Soft	8348	8350
Firm	8349	8351



sterile R skini CONTRA	144pk	500pk
Soft	8352	8354
Firm	8353	8355



sweep	144pk	500pk
Soft	8357	8359
Firm	8358	8360



TRUE INNOVATION  
FOR EVERY  
POLISHING NEED



It's a complete line for your practice. Our Ultrapro™ Tx Air handpiece is designed to be light, ergonomic, and comfortable. Developed to provide maximum comfort for you and your patient, it works with any DPA for high versatility. Our Ultrapro™ Tx Sweep DPA significantly reduces splatter—a truly innovative solution as part of the whole Ultrapro Tx product line.

# ultrapro™ Tx *air*

Designed with comfort in mind, both you and your patient will enjoy the smooth polishing experience. Lightweight and ergonomic, the Ultrapro Tx Air handpiece will be an efficient and effortless tool you'll love to use.

- Ergonomic shape for maximum comfort
- Lightweight (2.2 oz) to reduce hand and arm fatigue
- 360° swivel for easy rotation
- Soft-start motor to minimize splatter
- Minimal vibration for maximum comfort
- Cleans effectively at a maximum speed of 4,000 RPMs
- Universal E-type motor for multifunctionality
- Compatible with nearly all disposable prophylaxis angles



# ultrapro™ Tx *prophy paste*

Ultrapro™ Tx prophy paste and Ultrapro™ Tx Pure prophy paste perfectly round out all your cleaning needs. Low splatter and easy to rinse, they provide a comfortable cleaning experience. With a variety of flavors to choose from, Ultrapro Tx prophy pastes will impress your patients with their great taste and leave them looking forward to their next visit.



- Contains 1.23% fluoride ion
- Low splatter formula
- Available in Cool Mint, Walterberry™, Bubble Gum, and Orange Dreamsicle flavors
- Rinses easily and completely to eliminate residual grittiness
- Gluten free




- Free of fluoride, flavors, dyes, and oil
- Low splatter formula
- Rinses easily and completely to eliminate residual grittiness
- Ideal for use before bonding procedures
- Gluten free


# ultrapro™ Tx *disposable prophylaxis angles*

Ultrapro™ Tx disposable prophylaxis angles are innovative, effective, and comfortable to use. They feature smooth, quiet gears and an ergonomic design so that both you and your patient will have a comfortable experience. With an advanced cup design both inside and out, the Extra, Skini, and Sweep angles are built to clean better than ever. The Ultrapro Tx Sweep disposable prophylaxis angle features a new, innovative brush guard that prevents saliva and paste from collecting and roping around the cup to help reduce splatter up to 95%<sup>1</sup>.

- Innovative cup design for reduced splatter and efficient cleaning
- Outer ridges for improved interproximal cleaning (Extra, Skini, and Sweep angles)
- Designed for better access and improved visibility, with 20% shorter head and 25% slimmer neck (Skini and Sweep angles)
- Available in right-angle or contra-angle design (Extra and Skini angles)
- All designs feature optimal flare
- Ergonomic design for comfort
- Not made with natural rubber latex
- Available in soft or firm cup designs
- Skini angle is sterile to minimize risk of infection
- Sweep angle features an innovative brush design that significantly reduces splatter



1. Data on file.



*original*  
(Available in soft or firm cup)



*extra*  
(Available in soft or firm cup, right-angle or contra-angle)



*skini*  
(Available in soft or firm cup, right-angle or contra-angle)



*sweep*  
(Available in soft or firm cup)