

TOTAL CONTROL AT YOUR FINGERTIPS

That's what you'll get with Comfort Hub™,
the new design for our Metal Dento-Infusor™ tip.



Concave ergonomic hub for comfort

Greater wing surface area for easier control

The new Comfort Hub™ design has larger wings for easy attachment and removal
—providing a secure, gentle grasp of the tip.

Metal Dento-Infusor™ Tip

INTRODUCING THE NEW **COMFORT HUB™** DESIGN



- Attach to a syringe to apply any of Ultradent's tissue management products such as: **Astringedent™, Astringedent X™, ViscoStat™, and ViscoStat™ Clear hemostatics**
- The most recognized and utilized tip from Ultradent, and the first sold when the company was founded
- A unique, slightly concave ergonomic hub allows for easier attachment and removal
- The Comfort Hub design creates a greater surface area to make grip and torque lighter on the fingers



Metal Dento-Infusor™ Tip with Comfort Hub™

- Places hemostatic agents precisely and effectively removes superficial coagulum
- Blunt, bent cannula with padded brush enables gentle pressure in the sulcus
- Ultradent's first tip, the "MDI" remains paramount for successful tissue management

Designed for: Astringedent™, Astringedent™ X, ViscoStat™, ViscoStat™ Clear, and Peak™ Universal Bond.

LOK-TITE™ + COMFORT HUB™	100pk	500pk
19 ga Metal Dento-Infusor	4955	4956

