



**4670 - Jiffy Diamond Strips Perforated Narrow Assorted 15pk**  
5 x Each Strip Stainless Steel electroplated diamond strips  
(5 Medium, 5 Fine, and 5 Xfine)



**4674 - Jiffy Diamond Strips Perforated Wide Assorted 15pk**  
5 x Each Stainless Steel electroplated diamond strips  
(5 Medium, 5 Fine, and 5 Xfine)




**4671 - Jiffy Diamond Strips Perforated Narrow Medium 10pk**  
10 x Stainless Steel electroplated diamond strips



**4683 - Jiffy Diamond Strips Perforated Wide Medium 10pk**  
10 x Stainless Steel electroplated diamond strips



**4672 - Jiffy Diamond Strips Perforated Narrow Fine 10pk**  
10 x Stainless Steel electroplated diamond strips



**4676 - Jiffy Diamond Strips Perforated Wide Fine 10pk**  
10 x Stainless Steel electroplated diamond strips



**4673 - Jiffy Diamond Strips Perforated Narrow Xfine 10pk**  
10 x Stainless Steel electroplated diamond strips



**4677 - Jiffy Diamond Strips Perforated Wide Xfine 10pk**  
10 x Stainless Steel electroplated diamond strips

**4680 - Jiffy Proximal Saw Stainless Steel 10pk**  
10 x Stainless Steel ribbon saws



**DON'T SHORTCUT YOUR FINISH**

Ultradent's Jiffy™ composite finishing line provides a complete, quality solution to easily achieve a beautiful finish—without shortcut results.



SCAN TO VIEW PRICING

The Jiffy Natural composite finishing and polishing system consists of malleable spiral-shaped wheels that effectively polish all tooth surfaces.


- Easily re-creates the luster of natural enamel
- Complete procedural kit includes everything you need
- Ultradent diamond grit gives a beautiful finish on any composite material
- Ideal for finishing Ultradent's Transcend™ universal composite and Vit-l-escence™ restorative material
- Efficient 2-step polishing process saves you chair time



**6304-1 - Jiffy Natural Composite Polishing Kit**  
1 x Jiffy RA Medium Natural Composite 14 mm Wheel  
1 x Jiffy RA Fine Natural Composite 14 mm Wheel  
1 x Jiffy Medium Natural Composite Twirl  
1 x Jiffy Fine Natural Composite Twirl



**6089-1 - Jiffy Natural RA Medium 14 mm Polishing Wheel 3pk**  
**6090-1 - Jiffy Natural RA Fine 14 mm Polishing Wheel 3pk**



**6305-1 - Jiffy Natural Medium Twirl 3pk**  
**6306-1 - Jiffy Natural Fine Twirl 3pk**

The Jiffy Spin shaping and finishing disk system is a high-quality solution to shape and finish restorations, at a smart price, rapidly and efficiently.

- Ultra-thin disks are perfect for interproximal finishing and shaping. Every component of the Jiffy system will assist clinicians in providing the esthetic restoration patients want.
- Aluminum oxide coated disks range from coarse to fine for smoothing and polishing, PLUS something only the Jiffy Spin shaping and finishing disks can offer: Extra-Coarse diamond grit for rapid shaping and gross removal!

75pk	Extra-Coarse	Coarse	Medium	Fine
10 mm	5582	5572	5573	5574

75pk	Extra-Coarse	Coarse	Medium	Fine
14 mm	5583	5576	5578	5579

**Jiffy Spin Shaping & Finishing Disks**



**5570 - Jiffy Spin 10 mm Shaping & Finishing Disk Kit**

2 x Mandrels, 50 x Extra-Coarse Disks, 50 x Coarse Disks, 50 x Medium Disks, 50 x Fine Disks



**5571 - Jiffy Spin 14 mm Shaping & Finishing Disk Kit**

2 x Mandrels, 50 x Extra-Coarse Disks, 50 x Coarse Disks, 50 x Medium Disks, 50 x Fine Disks

**5581 - Jiffy Spin Mandrels 10pk**



The Jiffy Original composite finishing and polishing system is excellent for contouring, finishing, and polishing any composite material.

- Polishing cups feature a flared, flexible, thin wall design ideal for polishing near gingiva, labial surfaces, and easily covers all tooth surfaces
- Features unique polishing brushes that facilitate unsurpassed polishing
- Unique disk shape gives you leverage on working surfaces



**4254 - Jiffy Composite Adjusting & Polishing Kit**  
3 x Each cups, disks, and points  
(1 coarse, 1 medium, 1 fine)  
2 x Jiffy Brushes (1 regular, 1 pointed)

	Cups 12pk	Disks 12pk	Points 12pk
Coarse	7011	7015	7019
Medium	7012	7016	7020
Fine	7013	7017	7021

	Cups 10pk	Disks 10pk	Points 10pk
HiShine	7014	7018	7022



**7010 - Jiffy Composite Polishing Variety Pack**  
6 x Cups (2 coarse, 2 medium, 2 fine)  
3 x Disks (1 coarse, 1 medium, 1 fine)  
9 x Points (3 coarse, 3 medium, 3 fine)



**850 - Jiffy Regular Brush 10pk**  
**1009 - Jiffy Pointed Brush 10pk**



**1029 - Jiffy Goat Hair Brush 10pk**



**5540 - Diamond Polish Mint Syringe 0.5 µm 2pk**  
**5541 - Diamond Polish Mint Syringe 1 µm 2pk**  
1.2 ml syringes