

ORIGINAL



THE ONE YOU'VE LOVED FOR YEARS

The Omni-Matrix Original disposable retainer and matrix features a mechanism that locks the matrix into place before placement.

- 2201 Omni-Matrix Winged 0.001 Orange Refill
- 2202 Omni-Matrix Winged 0.0015 Green Refill
- 2204 Omni-Matrix Winged Pedo Purple Refill

48 x Omni-Matrix winged



- 2203 Omni-Matrix Winged Mylar Red Refill

48 x Omni-Matrix winged



- 1101 Omni-Matrix Wingless 0.001 Orange Refill
- 1102 Omni-Matrix Wingless 0.0015 Green Refill
- 1104 Omni-Matrix Wingless Pedo Purple Refill

48 x Omni-Matrix wingless



- 1103 Omni-Matrix Wingless Mylar Red Refill

48 x Omni-Matrix wingless



THE REDESIGNED
OMNI-MATRIX™
DISPOSABLE RETAINER AND MATRIX

YOUR CIRCUMFERENTIAL
BAND SOLUTION




SCAN
TO VIEW
PRICING

YOUR CIRCUMFERENTIAL BAND SOLUTION

The redesigned Omni-Matrix disposable retainer and matrix is a simple restorative tool designed to perfectly customize to any preparation. The band's circumference can be easily adjusted simply by twisting the handle and the pivoting head allows it to access any quadrant of the mouth. Once the restoration is complete, the Omni-Matrix band easily releases without disturbing the restorative material.

- Innovative shape allows procedural visibility and patient comfort
- Ultra-thin burnishable stainless steel adapts to all preparations
- Unique winged and wingless styles meet individual case needs
- Disposable design saves you time and money



WINGLESS

- 7701 Omni-Matrix Universal 6.5 mm Wingless
.001" (.025 mm) SS Orange
- 7702 Omni-Matrix Universal 6.5 mm Wingless
.0015" (.038 mm) SS Green
- 7704 Omni-Matrix Narrow 5.2 mm Wingless
.0015" (.038 mm) SS Purple

48 x Omni-Matrix wingless



- 7703 Omni-Matrix Universal 6.5 mm Wingless
.0025" (.064 mm) Mylar Red

48 x Omni-Matrix wingless



WINGED

- 8801 Omni-Matrix Universal 6.5 mm Winged
.001" (.025 mm) SS Orange
- 8802 Omni-Matrix Universal 6.5 mm Winged
.0015" (.038 mm) SS Green
- 8804 Omni-Matrix Narrow 5.2 mm Winged
.0015" (.038 mm) SS Purple

48 x Omni-Matrix winged



- 8803 Omni-Matrix Universal 6.5 mm Winged
.0025" (.064 mm) Mylar Red

48 x Omni-Matrix winged

