


ORIGINAL



THE ONE YOU'VE LOVED FOR YEARS

The Omni-Matrix Original disposable retainer and matrix features a mechanism that locks the matrix into position prior to placement.

WINGED

Stainless Steel	Winged 48pk
6.5 mm — .001" (.025 mm)	2201
6.5 mm — .0015" (.038 mm)	2202
5.2 mm — .0015" (.038 mm)	2204



Mylar	Winged 48pk
6.5 mm — .0025" (.064 mm)	2203



WINGLESS

Stainless Steel	Wingless 48pk
6.5 mm — .001" (.025 mm)	1101
6.5 mm — .0015" (.038 mm)	1102
5.2 mm — .0015" (.038 mm)	1104




Mylar	Wingless 48pk
6.5 mm — .0025" (.064 mm)	1103



THE REDESIGNED OMNI-MATRIX™ DISPOSABLE RETAINER AND MATRIX

YOUR CIRCUMFERENTIAL BAND SOLUTION






SCAN
TO VIEW
PRICING

YOUR CIRCUMFERENTIAL BAND SOLUTION

The redesigned Omni-Matrix disposable retainer and matrix is a simple restorative tool designed to perfectly customize to any preparation. The band's circumference can be easily adjusted simply by twisting the handle, and the pivoting head allows it to access any quadrant of the mouth. Once the restoration is complete the Omni-Matrix band easily releases without disturbing the restorative material.

- Innovative shape allows procedural visibility and patient comfort¹
- Ultra-thin burnishable stainless steel adapts to all preparations¹
- Unique winged and wingless styles meet individual case needs
- Disposable design saves you time and money



WINGLESS

Stainless Steel	Wingless 48pk
6.5 mm — .001" (.025 mm)	7701
6.5 mm — .0015" (.038 mm)	7702
5.2 mm — .0015" (.038 mm)	7704



Mylar	Wingless 48pk
6.5 mm — .0025" (.064 mm)	7703



WINGED

Stainless Steel	Winged 48pk
6.5 mm — .001" (.025 mm)	8801
6.5 mm — .0015" (.038 mm)	8802
5.2 mm — .0015" (.038 mm)	8804



Mylar	Winged 48pk
6.5 mm — .0025" (.064 mm)	8803

