


CUSTOMIZE PERFECT PREPARATIONS YOUR WAY:

THE ULTIMATE TRIODENT BUYING GUIDE



SuperCurve Non-Stick Matrices

- Ultra-thin non-stick finish minimizes adhesion to the bonding agent or composite to the matrix band
- Matrix shape features anatomical accuracy and curvature
- Gingival apron provides matrix coverage for deeper preparations



V3 Tab-Matrices

- Matrix bands feature significantly greater horizontal curvature ensuring that more curvature remains as the matrix wraps around the tooth
- Gingival apron provides matrix coverage for deeper preparations



ClearMetal™ Matrix

- Resin-laminated micro-windows allow light to pass through the matrix
- Greater matrix curvature hugs tooth for stability during placement
- Non-stick finish makes matrix removal easy



Tab-Matrices

- Matrix bands are super thin
- Gingival apron provides matrix coverage for deep preparations

CUSTOMIZE PERFECT PREPARATIONS YOUR WAY

THE ULTIMATE TRIODENT BUYING GUIDE




Wave-Wedges™

- Wave-shaped wings compress on entry and flare again upon exit for easy placement and a true fit and seal
- The hollow underside allows another wedge to be placed securely from the opposite side
- Wedges are designed to minimize impingement of soft tissue




Adaptive Wedge™

- All of the features of the Wave-Wedge™ with the added benefit of transparency, making it possible for you to initiate curing through the wedge
- Notches enhance flexibility for precise matrix adaptation and retention



V-Wedges

- V-Wedges are rigid enough to provide tooth separation
- The wedges are self-guided and stackable, which make them easy to place
- The wedges also adapt to the contours of the interproximal space



V3 Rings

Universal and Narrow

- Super-elastic NiTi gives outstanding strength and consistent tooth separation force
- V-Shape tines accommodates the wedge from both sides and prevents the ring from collapsing into wide preparations
- Stackable rings for multiple simultaneous restorations