



Umbrella	5pk	20pk	40pk
Soft Tissue Retractor Medium	4870	4871	5162



**352 - OraSeal Syringe Kit**  
2 x 1.2 ml OraSeal Caulking syringes  
2 x 1.2 ml OraSeal Putty syringes  
4 x Black Mini tips  
20 x White Mac tips



1.2 ml syringe	4pk	20pk
Caulking	351	354
Putty	353	355



**4100 - PropGard Kit**  
5 x Each PropGards large yellow and regular purple  
10 x Each Tongue guards left and right



PropGard	Regular 5pk	Large 5pk
Soft	4105	4103
Firm	4101	4104



**4102 - Tongue Guard 10pk**  
5 x Each right and left



**311 - DermaDam Medium 0.20 mm 36pk**  
6"x 6" (15 cm x 15 cm)



**314 - DermaDam Heavy 0.25 mm 36pk**  
6"x 6" (15 cm x 15 cm)



**299 - DermaDam Medium Synthetic 0.20 mm 20pk**  
6"x 6" (15 cm x 15 cm)

Reliable Moisture Control,  
Clear Access, Superior Results

Keep Patients  
Comfortable



## Umbrella™ Soft Tissue Retractor

The new Umbrella retractor was developed to make patient comfort its key component:

- 💧 A new, innovative tongue-retraction design allows the tongue to comfortably rest behind the tongue guard
- 💧 Easy to place
- 💧 Naturally and gently keeps the mouth propped open, which means no pulling or stretching of the lips
- 💧 Designed with anatomically placed/shaped bumpers, so clinicians can rest a hand on the patient's mouth without any pain points
- 💧 Provides relief to gaggers—for many it doesn't initiate the gag reflex



## PropGard™ Mouth Prop

Our PropGard mouth prop protects the tongue from trauma and prevents adjacent teeth, tongue, and mouth closure contamination.

- 💧 Achieve a dry working environment



## OraSeal™ Caulking and Putty

OraSeal™ Caulking and Putty materials effectively adhere to a wet dental dam, tissue, teeth, and metal. Even under water or saliva, OraSeal materials excel at preventing leakage during treatment, protect the gingiva from exposure to peroxide or hydrofluoric acid, and block out unwanted spaces for impressions.

- 💧 In conjunction with the tear-resistant DermaDam rubber or DermaDam™ Synthetic dental dam, it's easy to achieve a dry working environment.

## DermaDam™ Rubber Dam

DermaDam rubber dental dam is made from powder-free pure latex rubber. It has a low dermatitis potential and the lack of powder reduces the occurrence of allergic reactions.

- 💧 DermaDam Synthetic dental dam is not made with natural rubber latex but is designed to be just as flexible and durable.

