

# Perfect Delivery. Optimal Performance. Easy Removal.



simplifying



endodontics

## UltraCal<sup>®</sup> XS and Citric Acid 20%

**UltraCal XS** calcium hydroxide paste (pH 12.5) is easily delivered with the NaviTip<sup>®</sup> tip exactly where it is needed in the canal. Its aqueous consistency and radiopaque properties make it ideal as an endodontic treatment dressing for two-appointment RCTs.

When it comes time to remove UltraCal XS paste from the canal, look no further than **Ultradent<sup>®</sup> Citric Acid 20% solution**, delivered with the NaviTip<sup>®</sup> FX<sup>®</sup> tip. Citric Acid 20% solution dissolves calcium hydroxide, and the small fibers attached to the NaviTip FX tip easily scrub the walls of the canal, removing any remaining calcium hydroxide as well as the smear layer.



**ULTRADENT**  
PRODUCTS, INC.

800.552.5512 [ULTRADENT.COM](http://ULTRADENT.COM)



## UltraCal<sup>®</sup> XS

30%–35% CALCIUM HYDROXIDE PASTE

Certified  
**GF**  
Gluten-Free



## NaviTip<sup>®</sup> Tips



20pks	29 ga	30 ga
17 mm - White	5112	1249
21 mm - Yellow	5113	1349
25 mm - Blue	5114	1250
27 mm - Green	5115	1354
Assorted	5116	1351

50pks	29 ga	30 ga
17 mm - White	1378	1421
21 mm - Yellow	1374	1422
25 mm - Blue	1376	1423
27 mm - Green	1377	1424
Assorted	1379	3319



**1027 UltraCal XS Kit**  
4 x 1.2 ml syringes  
20 x NaviTip tips  
(Assorted 29 ga - 5 of each length)



**606 UltraCal XS Refill**  
4 x 1.2 ml syringes



**5117 UltraCal XS Econo Refill**  
20 x 1.2 ml syringes



**329 Citric Acid 20% IndiSpense Syringe**  
1 x 30 ml syringe



## NaviTip<sup>®</sup> FX<sup>®</sup> Tip

1452 17 mm - 20pk



## NaviTip<sup>®</sup> FX<sup>®</sup> Tip

1454 25 mm - 20pk

