

NO NEED TO RUSH

You have 5 full days.



Chromaclone™
5-day stability alginate

Chromaclone™

Chromaclone™ alginate

The Freedom of Five Full Days

In your busy dental office, every second counts. Every extra minute you save lets you breathe a little easier. With Chromaclone 5-day alginate impression material you'll get more than a few extra minutes—you'll get a few extra days. That's because Chromaclone alginate's formulation offers a full five days to pour your impressions, not just the few hours you're used to.*

And as always, Chromaclone alginate offers the same dust-free formula and reliable set times you depend on. The Fast set time (53 seconds) features a convenient color-changing formula, while the Extra Fast set time (45 seconds) is ideal for more sensitive patients. And both are available in a delicious Bubble Gum flavor. **YUM!**



Fast Set: 53 seconds time in mouth

4033 - Chromaclone 5-day FAST SET Kit
1 x 1 lb bag of Chromaclone 5-day
1 x Plastic container, water cup, and scoop



4034 - Chromaclone 5-day FAST SET Refill
1 x 1 lb bag of Chromaclone 5-day



Extra Fast Set: 45 seconds time in mouth

5502 - Chromaclone Ortho 5-day EXTRA FAST SET Kit
1 x 1 lb bag of Chromaclone 5-day
1 x Plastic container, water cup, and scoop



5503 - Chromaclone Ortho 5-day EXTRA FAST SET Refill
1 x 1 lb bag of Chromaclone 5-day



*Wrap the impression with a damp paper towel and store in a sealed plastic bag at room temperature (23°C/73°F) to ensure dimensional stability for up to five days.



SCAN
TO VIEW
PRICING