

# SAFEGUARD

your patients with **THE** orthodontic primer and sealant in one



**opal seal**<sup>™</sup>   
the ortho primer/sealant

*The revolutionary 38% filled Opal Seal primer and sealant releases and recharges fluoride.*

## RELEASING AND RECHARGING FLUORIDE

### Prime

- Drying agent chases moisture and draws fluoride-containing resin into etched enamel
- Single-step application ensures stronger bonds
- 38% filled with substantial glass ionomer and nanofillers for penetration of fissures to offer strong, long-lasting coverage and mechanical retention
- Highly filled layer
- Nonyellowing formulation cures translucent and resists staining

### Power Up

- Releases fluoride and recharges uptake
- Fluorescing properties are detectable under a UV black light to facilitate removal and reapplication
- Compatible with any bonding system



Opal Seal sealant is translucent after curing and resists staining.



Opal Seal sealant's fluorescing properties make it easy to remove under a UV black light.

### METAL OPAL™ SEAL TIP

The soft fibers and metal cannula are ideal for applying a thin, even layer to etched enamel.



### 500-060 - Opal Seal Syringe Intro Kit

- 1 x 1.7 g Opal Seal syringe
- 1 x 1.5 g Opal Etch syringe
- 10 x Metal Opal Seal tips
- 10 x Opal Blue tips



### 500-063 - Opal Seal Syringe Sample Kit

- 1 x 0.6 g Opal Seal syringe
- 3 x Metal Opal Seal tips



### 500-040 - Opal Seal Adhesives Syringe Kit

- 1 x 1.7 g Opal Seal syringe
- 1 x 0.9 g Opal Bond MV syringe
- 1 x 2.3 g Opal Bond Flow syringe
- 1 x 1.5 g Opal Etch syringe
- 20 x Metal Opal Seal tips
- 20 x Opal Blue tips
- 20 x Opal Bond Flow tips



### 500-061 - Opal Seal Syringe Kit

- 4 x 1.7 g Opal Seal syringes
- 60 x Metal Opal Seal tips



### 500-062 - Opal Seal Syringe 20pk

- 1.7 g syringes

