

VALO™

ORTHOD

CORDLESS LED CURING LIGHT




ESPECIALLY FOR ORTHODONTICS



The VALO™ Ortho Cordless broadband LED curing light was designed specifically for orthodontists.


- Unique Xtra Power Quadrant mode quickly and efficiently cures five teeth with one touch of a button
- Powerful, optimally collimated beam delivers a consistent, even cure directly over the labial face of the bracket
- Slim, low-profile design allows unprecedented access anywhere in the mouth
- Intuitive timer controls are easy to operate and user-friendly
- Highly efficient LEDs keep the wand body cool to the touch
- Battery-operated, cordless wand provides optimal convenience and flexibility
- Rechargeable batteries are environmentally responsible, safe, and inexpensive
- Comes with an exclusive 5-year manufacturer's warranty



The Black Light lens is used to detect fluorescent particles in resins for easy differentiation from natural enamel.



The Black Light lens magnetically attaches to the VALO Ortho light.



Ultra-high-energy broadband LEDs produce a high-intensity light at 385–515 nm.

5942 - VALO Ortho Cordless Kit

- 1 x VALO Ortho Cordless LED curing light
- 4 x Rechargeable batteries
- 1 x Battery charger
- 1 x Charging unit power supply
- 1 x Handpiece mounting bracket holder
- 1 x VALO cordless light shield
- 1 x Black Light lens
- 1 x Sample pack barrier sleeves

ACCESSORIES:

- 5962 - Battery Charger
- 5963 - Rechargeable Batteries (2 qty)
- 5961 - Power Supply with Universal Plugs for Battery Charger
- 5939 - Black Light Lens (1 qty)
- 5929 - VALO Cordless Light Shield
- 1667 - Surface Mounting Bracket
- 508 - UltraTect™ Orange-Lens Glasses

- 4667 - Barrier Sleeve Refill (100 qty)



SCAN
TO VIEW
PRICING

