

Unlimited Movement

 **ultrapro™ Tx**
cordless

Meet the new Ultrapro™ Tx Cordless handpiece—the hygiene handpiece that gives you the power of a corded handpiece with the unlimited movement of cordless. The cordless design means you have the ability to approach a tooth from any angle with no cord drag.



Prophylaxis Made Easy

- Cordless design allows for unlimited movement
- Compatible with most DPAs for maximum usability
- A single wireless charge can last an entire day
- Automatic shutoff after five minutes
- Working output of up to 4,000 RPM
- Wireless Bluetooth foot control
- Lightweight and ergonomic
- 2-year warranty



5275 - Ultrapro Tx Cordless Propphy Motor **Kit**

- 1 x Cordless propphy motor
- 3 x Sheaths
- 1 x Foot control
- 1 x Charger base
- 2 x Wall chargers
- 4 x Propphy angles
- 100 x Barrier sleeves
- 6 x Black sheath O-rings
- 2 x Green motor O-rings

5532 - Ultrapro Tx Cordless Propphy O-Ring

Replacement Part

- 6 x Black sheath O-rings
- 2 x Green motor O-rings



5533 - Ultrapro Tx Cordless Propphy Foot Control

Replacement Part

- 1 x Foot control



5276 - Ultrapro Tx Cordless Propphy Sheath **Kit**

- 1 x Sheath
- 6 x Black sheath O-rings
- 2 x Green motor O-rings



5534 - Ultrapro Tx Cordless Propphy Charger

Base Replacement Part

- 1 x Charger base



\$119.00 ea | \$113.05 ea 3/more | \$107.10 ea 5/more

5531 - Ultrapro Tx Cordless Propphy Wall Charger

Replacement Part

- 1 x Wall Charger



5499 - Ultrapro Tx Cordless Propphy Motor Barrier

Sleeves 500pk

- 500 x Barrier sleeves



SCAN
TO VIEW
PRICING

