

4243 Jiffy Composite Finishing Kit
 3 x Jiffy T&F Flames (1 7901, 1 Long, 1 Long 30-Blade)
 2 x Jiffy Cylinders (1 T&F Beveled, 1 7572 Flat)
 2 x Jiffy T&F Eggs (1 7406 Small, 1 7408 Large)
 1 x Jiffy Interproximal Diamond



4245 Jiffy 7901 T&F Flame 5pk
4251 Jiffy T&F Long Flame 5pk
4252 Jiffy T&F Long Flame 30-blade 5pk
 5 x Flames



4244 Jiffy 7572 Flat Cylinder 5pk
4249 Jiffy T&F Beveled Cylinder 5pk
 5 x Cylinders



4246 Jiffy 7408 T&F Egg Large 5pk
4247 Jiffy 7406 T&F Egg Small 5pk
 5 x Eggs



4248 Jiffy Interproximal Diamond 5pk
 5 x Interproximal Diamonds



Contact your Ultradent sales representative to find out more about Jiffy® finishing systems.



Jiffy
 COMPOSITE SYSTEM



NEW
Jiffy natural
 UNIVERSAL
 COMPOSITE
 SYSTEMS

DON'T SHORTCUT YOUR FINISH

Ultradent's Jiffy® Composite Finishing Line provides a complete, quality solution to easily achieve a beautiful finish—without shortcut results.



NEW

The new Jiffy Natural Universal Composite finishing and polishing systems consist of malleable spiral-shaped wheels that naturally conform to all tooth surfaces.

- Easily re-creates the luster of natural enamel
- Complete procedural kit includes everything you need
- Ultradent diamond grit gives a beautiful finish on any composite material
- Ideal for finishing Ultradent's Mosaic[®] universal composite

6104 **NEW** Jiffy Natural Composite Carbide Finishing & Polishing Kit

- 1 x Jiffy RA Medium Natural Composite 14 mm Wheel
- 1 x Jiffy RA Fine Natural Composite 14 mm Wheel
- 1 x Jiffy RA Medium Composite Cup
- 1 x Jiffy RA Fine Composite Cup
- 1 x Jiffy TF6 8 BLD Needle
- 1 x Jiffy 7901 T&F Flame
- 1 x Jiffy 7408 T&F Egg Large
- 1 x Jiffy Interproximal Diamond



6089 Jiffy RA Medium Natural Composite 14 mm Wheel 3pk
6090 **NEW** Jiffy RA Fine Natural Composite 14 mm Wheel 3pk
3 x Wheels



6109 **NEW** Jiffy RA Assorted Natural Composite 14 mm Wheel 6pk
6 x Wheels



6098 **NEW** Jiffy TF6 8 BLD Needle 5pk
5 x Needles



4248 Jiffy Interproximal Diamond 5pk
5 x Interproximal Diamonds



6105 **NEW** Jiffy Natural Composite Diamond Finishing & Polishing Kit

- 1 x Jiffy RA Medium Natural Composite 14 mm Wheel
- 1 x Jiffy RA Fine Natural Composite 14 mm Wheel
- 1 x Jiffy RA Medium Composite Cup
- 1 x Jiffy RA Fine Composite Cup
- 1 x Jiffy D6F Fine Needle
- 1 x Jiffy D8F Fine Flame
- 1 x Jiffy Diamond Fine Egg Large
- 1 x Jiffy Interproximal Diamond



6087 Jiffy D8F Fine Flame 5pk
6088 Jiffy Diamond Fine Egg Large 5pk
6099 **NEW** Jiffy D6F Fine Needle 5pk
5 x Each



4245 Jiffy 7901 T&F Flame 5pk
4246 Jiffy 7408 T&F Egg Large 5pk
5 x Each



838 Jiffy Medium Cup 20pk
841 Jiffy Fine Cup 20pk
20 x Cups



The original Jiffy Composite finishing and polishing system is excellent for contouring, finishing, and polishing any composite material.

- Polishing cups feature a flared, flexible thin-wall design ideal for polishing near the gingiva
- Features unique polishing brushes that facilitate unsurpassed polishing
- Ultradent formulation gives a beautiful finish on any composite material

4254 Jiffy Composite Adjusting & Polishing Kit

- 2 x Jiffy Brushes (1 Regular, 1 Pointed)
- 3 x Jiffy Cups (1 Coarse, 1 Medium, 1 Fine)
- 3 x Jiffy Disks (1 Coarse, 1 Medium, 1 Fine)
- 3 x Jiffy Points (1 Fine each)



890 Jiffy Coarse Cup 20pk
838 Jiffy Medium Cup 20pk
841 Jiffy Fine Cup 20pk
20 x Cups



891 Jiffy Coarse Disk 20pk
840 Jiffy Medium Disk 20pk
843 Jiffy Fine Disk 20pk
20 x Disks



892 Jiffy Coarse Point 20pk
839 Jiffy Medium Point 20pk
842 Jiffy Fine Point 20pk
20 x Points



848 Jiffy Polishing Variety Pack

- 5 x Cups (2 Coarse, 2 Medium, 1 Fine)
- 5 x Disks (2 Coarse, 2 Medium, 1 Fine)
- 10 x Points (4 Coarse, 4 Medium, 2 Fine)



3060 Jiffy HiShine Points 10pk
3061 Jiffy HiShine Cups 10pk
3062 Jiffy HiShine Disks 10pk
10 x HiShines



850 Jiffy Regular Brush 10pk
1009 Jiffy Pointed Brush 10pk
10 x Brushes



1029 Jiffy Goat Hair Brush 10pk
10 x Brushes



5540 Ultradent Diamond Polish Mint 0.5 µm Mini Refill
5541 Ultradent Diamond Polish Mint 1 µm Mini Refill
2 x 1.2 ml syringes



Systems are available with autoclavable aluminum blocks. Ultradent recommends the use of a block when autoclaving.