

speedhand solution for below & above elbow. Order form.

Customer

Customer no.	
Company	
Street	
Postal code	City
O&P professional	

Shipping address (if different from customer address)

Customer no.	
Company	
Street	
Postal code	City
Commision	

Preferred shipping date _____

(Generally, delivery takes place 1 – 2 days after our dispatch)

Prosthesis user

Name _____

Order-ID _____

amputation level transradial (below elbow)
transhumeral or higher (above elbow)

Fitting

Initial prosthetic fitting

User was not previously fitted with a myoelectric hand.

Follow-up prosthetic fitting

User was previously fitted with a myoelectric hand

Serial-no. previous hand _____

amputation side left right

Terminal device *speedhand*

Variations *speedhand*

Art. no.	Side	Size	thumb sensor	
8E2=L-XS	left	XS (7)	0 - without	1 - with
8E2=L-S	left	S (7 ¼)	0 - without	1 - with
8E2=L-M	left	M (7 ¾)	0 - without	1 - with
8E2=L-L	left	L (8 ¼)	0 - without	1 - with
8E2=R-XS	right	XS (7)	0 - without	1 - with
8E2=R-S	right	S (7 ¼)	0 - without	1 - with
8E2=R-M	right	M (7 ¾)	0 - without	1 - with
8E2=R-L	right	L (8 ¼)	0 - without	1 - with

Warranty packages *speedhand*

Art. no.	Description
SP-8E2=3	3 years warranty package
SP-8E2=3+3	3 years warranty package extension
SP-8E2=6	6 years warranty package

Wrist connection (modular wrists)

Art. no.	Description
10V61=1	wrist flex
10V62=1	wrist eqd
10V63=1	wrist short
10V64=1	wrist transcarpal

Glove covers

Variations	Gender
Skin Natural (-N)	female
Standard	male

Glove color

Skin Natural (-N)									
Color	2	4	6	8	11	14	16		
Standard prosthetic glove									
Color	1	2	3	4	5	6	7	8	9
	10	11	12	13	14	15	16	17	18

Control unit

Art. no.	Description
13E522=1	myosmart

connectgrip app



speedhand solution for below & above elbow. Order form.

Terminal device *DynamicArm*

Fitting

Initial prosthetic fitting

User was not previously fitted with a *DynamicArm/DynamicArm plus*.

Follow-up prosthetic fitting

User was previously fitted with a *DynamicArm/DynamicArm plus*.

Serial-no. _____

Warranty packages *DynamicArm*

Art. no.	Description
12K100N=3	3 years warranty package
12K100N=5	5 years warranty package
12K100N=3+2	3 years warranty package extension

Electrodes

Art. no.	Description	Filter	Qty.		
13E202=50	Suction socket electrode	50 Hz	2	3	4
13E202=60	Suction socket electrode	60 Hz	2	3	4
13E200=50	Standard electrode	50 Hz	2	3	4
13E200=60	Standard electrode	60 Hz	2	3	4
13E401 + 13E400=50	Remote electrodes*	50 Hz	2	3	4
13E401 + 13E400=60	Remote electrodes*	60 Hz	2	3	4

* 2x 13E401 are configured automatically. Remote electrodes can only be used in combination with 13E522 myosmart.

Electrode cable

Art. no.	Cable length	Qty.		
13E129=G300	300 mm	2	3	4
13E129=G600	600 mm	2	3	4
13E129=G1000	1000 mm	2	3	4
13E401/13E402=90	90 mm	2	3	4
13E401/13E402=140	140 mm	2	3	4

Note: for each 13E400=* 2 domes and 13E401=* 3 domes are required.

Electrode domes for 13E400/401

Art. no.	Height	Qty.
13Z161	low	
13Z162	medium	
13Z163	high	

ErgoArm Hybrid Plus

Art. no.	Diameter	Color	
12K44=45	45 mm	(4) skin	(7) black
12K44=50	50 mm	(4) skin	(7) black
		(11) brown*	(15) dark brown*

ErgoArm Electronic Plus

Art. no.	Diameter	Color	
12K50=45	45 mm	(4) skin	(7) black
12K50=50	50 mm	(4) skin	(7) black
		(11) brown*	(15) dark brown*

DynamicArm

Art. no.	Diameter	Color	
12K100N=45	45 mm	(4) skin	(7) black
12K100N=50	50 mm	(4) skin	(7) black
		(11) brown*	(15) dark brown*

* Only for elbows with diameter 50mm

Rotation options with EQD connection

Description	Consists of
Passive rotation	Coupling piece 10S4 and Coaxial plug 9E169
Active rotation	Electric wrist rotator 10S17 and MyoRotronic 13E205
Active rotation	Electric wrist rotator 10S17, when myosmart or DynamicArm is selected

Note: MyoRotronic is not shipped if myosmart is selected, except in combination with ErgoArm and „Active rotation“

Lamination ring

Art. no.	Diameter		
10S1	45 mm	50 mm	54 mm

Energy-management and charger unit

Art. no. battery	Battery design	Art. no. w. charger	Color
757B35=3	MyoEnergy Integral	757L35	–
757B35=4	MyoEnergy Integral	757L35	–
757B35=5	MyoEnergy Integral	757L35	–
757B20	Energy Pack	757L20	(4) skin
757B20-1	Energy Pack	757L20	(11) brown
757B20-2	Energy Pack	757L20	(15) dark brown
757B20-3	Energy Pack	757L20	(7) black
757B21	Energy Pack	757L21	(4) skin
757B21-1	Energy Pack	757L21	(11) brown
757B21-2	Energy Pack	757L21	(15) dark brown
757B21-3	Energy Pack	757L21	(7) black

Note: The charger will be selected automatically.

speedhand solution for below elbow. Order form.

Optional accessories

Skeo up daily assist

Art. no.	Size					Qty.
14Y4	140	160	180	200	220	

Greifer

Art. no.	Qty.
8E33=9-1 DMC VariPlus	
8E34=9-1 DMC VariPlus	

Accessories socket

Art. no.	Description
12V10	Suction socket valve
99B13=21	PVC-pull-in tube
99B13=21-7	PVC-pull-in tube, black
OC1560=TR	Easy Fit Arm for transradial fittings

Tool set for wrist assembly

Art. no.	Description
711M106=1	Mounting plate for wrist assembly
709H11=1	Spanner wrench for wrist assembly

Optional accessories elbow

Textile harness

Art. no.	Side	Size	Qty.
21A47=L-S	left	S	
21A47=L-M	left	M	
21A47=R-S	left	L	
21A47=R-M	right	S	
21A47=R-L	right	M	
21A47=L-L	right	L	
OC1560=TH	Easy Fit Arm for transhumeral fittings.		

Accessoires myosmart

Art. no.	Description
13E101=1	Optional Beeper for myosmart control
13E102=1	Elbow adapter cable for DynamicArm and ErgoArm 12K50

Optional accessories for gloves

Art. no.	Description
640F12	Special cleaner
640F13	Pump sprayer for 640F12 cleaner

USB-Adapter for charging (additional, for travelling)

757L43

Note: The 757L35 charger is not included in the scope of delivery of the 757L43 and needs to be ordered separately for a complete second charger with USB port must be ordered separately.

Myo evaluation tool

757M20-2=50 myosmart cuff
757M20-2=60 myosmart cuff

Accessoires DynamicArm

Art. no.	Description
646C42=V1.6	Elbowsoft 1.6
60X5	BionicLink
13E100	AnalogAdater for DynamicArm

Carbon Cover

Art. no.	Description
13Z157=45	for DynamicArm (diameter 45mm)
13Z157=50	for DynamicArm (diameter 50mm)

Date _____ Place _____ Signature _____