

Maverick Xtreme AT.

Your strong partner.

The **Maverick Xtreme AT** was designed for highly active users who navigate varied indoor and outdoor environments and place a high value on effortless walking and the ability to go wherever life takes them.

Designed for optimal outcomes

- The fibreglass design delivers a high degree of durability and flexibility without limiting energy return
- Developed for high loads
- Split forefoot and heel spring for even better terrain adaptation
- A full length, unbolted sole plate ensures a smooth stance phase rollover
- Protected against fresh, salt and chlorinated water

More than a foot. A foundation.



Technical data

Mobility grade / K-level	3, 4
Max. body weight	166 kg / 365 lbs
Side	left (L), right (R)
Sizes	23 – 30 cm
Weight without footshell*	520 g / 18.34 oz
Footshell colour	beige, light brown
Weight with footshell*	820 g / 28.92 oz
Heel height	10 mm / 3/8"
Build height*	170 mm / 6 11/6"

* Reference size is 26 cm.

ottobock.



Smooth rollover
Robust and durable
even under
heavy loads
Supports dynamic
movements

F21

Order form.

Customer

Customer no.	
Company	
Street	
Postal code	City
O&P professional	

Shipping address (if different from customer address)

Customer no.	
Company	
Street	
Postal code	City
Customer ref.	

Selection of the spring stiffness relative to

1. body weight and activity as well as

Body weight (kg / lbs)	Activity level		
	Low	Moderat	High
44 – 52 / 100 – 115	1	1	2
53 – 59 / 116 – 130	1	2	3
60 – 68 / 131 – 150	2	3	4
69 – 77 / 151 – 170	3	4	5
78 – 88 / 171 – 195	4	5	6
89 – 100 / 196 – 220	5	6	7
101 – 116 / 221 – 255	6	7	8
117 – 130 / 256 – 285	7	8	9
131 – 147 / 286 – 325	8	9	–
148 – 166 / 326 – 365	9	–	–

2. foot size

Stiffness	1	2	3	4	5	6	7	8	9
Foot size (cm)									
23							–	–	–
24								–	–
25								–	–
26	–								–
27	–	–							
28	–	–						**	**
29	–	–	–	**	**	**	**	**	**
30	–	–	–	**	**	**	**	**	**

** Do not combine this configuration with a 3B5-4 Genium X4 for body weight > 125 kg (275 lbs)

Example for ordering: F21=L26-4-P/4

Contains: F21-00-04A26-SL Foot module
FTC-2A-126L4-SL Beige footshell

Quantity	Article no.	Side	Size	Stiffness	Colour
	F21	=	-	-	P /
		Side	Size (cm)	Stiffness	Colour
		R Right	23, 24, ..., 30	1, 2, ..., 9	4 Beige
		L Left			15 Light brown

Alternative ordering option – Article number foot module (without footshell)

Quantity	Article no.	Stiffness	Size	Side
	F21 - 00 - 0	A	-	S
		Stiffness	Size (cm)	Side
		1, 2, ..., 9	23, 24, ..., 30	R Right
				L Left

Alternative ordering option – Article number footshell

Quantity	Article no.	Size	Colour	Side
	FTC - 2 A - 1		4	S
		Size (cm)	Colour	Side
		23, 24, ..., 30	L Beige	R Right
			D Light brown	L Left



The scope of delivery includes the F21* foot module, the FTC-2A-1* footshell, the 2F11=1 heel wedge set and a black Spectra sock. Additionally, a tool for changing the footshell can be ordered (ACC-00-10300-00).

Compatible mechatronic knee joint



Genium X4