

micelangelo and *myosmart plus* for below & above elbow. Order form.

Customer

Customer no.	
Company	
Street	
Postal code	City
O&P professional	

Shipping address (if different from customer address)

Customer no.	
Company	
Street	
Postal code	City
Commision	

Preferred shipping date

(Generally, delivery takes place 1–2 days after our dispatch)

Prosthesis user

Name

Order-ID

amputation level ☐ transradial (below elbow)
☐ transhumeral or higher (above elbow)

Fitting

☐ **Initial prosthetic fitting**
 User was not previously fitted with a myoelectric hand.

☐ **Follow-up prosthetic fitting**
 User was previously fitted with a myoelectric hand.

Serial-no. previous hand

amputation side ☐ left ☐ right

Terminal device *micelangelo*

Variations *micelangelo*

Art. no.	Side	Size
8E5=L-S	left	S
8E5=L-M	left	M
8E5=R-S	right	S
8E5=R-M	right	M

Warranty packages *micelangelo*

Art. no.	Description
☐ SP-8E5=3	3 years warranty package
☐ SP-8E5=3+3	3 years warranty package extension (3 to 6 years)
☐ SP-8E5=6	6 years warranty package

Wrist connection (modular wrists)

Art. no.	Description
☐ 10V61=1	wrist flex
☐ 10V62=1	wrist eqd
☐ 10V63=1	wrist short
☐ 10V64=1	wrist transcarpal

Glove covers PVC

Gender
☐ female ☐ male

Glove color
☐ Skin Natural
 Color ☐ 2 ☐ 4 ☐ 6 ☐ 8 ☐ 11 ☐ 14 ☐ 16
☐ black ☐ transparent

Glove covers Silicone (only male and for handsizes M)

Glove color
☐ Color ☐ 2 ☐ 4 ☐ 6 ☐ 8 ☐ 10 ☐ 12 ☐ 14

Control unit

Art. no.	Description
☐ 13E522=1	myosmart
☐ 13E523=TR-1	myosmart plus TR
☐ 13E523=TH-1	myosmart plus TH

connectgrip app



micelangelo and myosmart plus **for below & above elbow.** **Order form.**

Elbow component *DynamicArm*

Fitting

☐ **Initial prosthetic fitting**
User was not previously fitted with a *DynamicArm*.

☐ **Follow-up prosthetic fitting**
User was previously fitted with a *DynamicArm*.

Serial-no. _____

Warranty packages *DynamicArm*

Art. no.	Description
☐ SP-12K100N=3	3 years warranty package
☐ SP-12K100N=5	5 years warranty package
☐ SP-12K100N=3+2	3 years warranty package extension (3 to 5 years)

Electrodes

Art. no.	Description	Filter	Qty.
☐ 13E202=50	Suction socket electrode	50 Hz	☐ 2 ☐ 3 ☐ 4 ☐ 5 ☐ 6 ☐ 7 ☐ 8
☐ 13E202=60	Suction socket electrode	60 Hz	☐ 2 ☐ 3 ☐ 4 ☐ 5 ☐ 6 ☐ 7 ☐ 8
☐ 13E200=50	Standard electrode	50 Hz	☐ 2 ☐ 3 ☐ 4 ☐ 5 ☐ 6 ☐ 7 ☐ 8
☐ 13E200=60	Standard electrode	60 Hz	☐ 2 ☐ 3 ☐ 4 ☐ 5 ☐ 6 ☐ 7 ☐ 8
☐ 13E400/401=*-50	Remote electrodes*	50 Hz	☐ 2 ☐ 3 ☐ 4 ☐ 5 ☐ 6 ☐ 7 ☐ 8
☐ 13E400/401=*-60	Remote electrodes*	60 Hz	☐ 2 ☐ 3 ☐ 4 ☐ 5 ☐ 6 ☐ 7 ☐ 8

* 2x 13E401 are configured automatically. Remote electrodes can only be used in combination with 13E522, 13E523=*

* number of electrodes: myosmart 2-4, myosmart plus 4-8

* mixing 13E20* electrodes and 13E40* electrodes is not possible

Electrode cable

Art. no.	Cable length	Qty.
☐ 13E129=G300	300 mm	☐ 2 ☐ 3 ☐ 4 ☐ 5 ☐ 6 ☐ 7 ☐ 8
☐ 13E129=G600	600 mm	☐ 2 ☐ 3 ☐ 4 ☐ 5 ☐ 6 ☐ 7 ☐ 8
☐ 13E129=G1000	1000 mm	☐ 2 ☐ 3 ☐ 4 ☐ 5 ☐ 6 ☐ 7 ☐ 8
☐ 13E400/401=G90-*	90 mm	☐ 2 ☐ 3 ☐ 4 ☐ 5 ☐ 6 ☐ 7 ☐ 8
☐ 13E400/401=G140-*	140 mm	☐ 2 ☐ 3 ☐ 4 ☐ 5 ☐ 6 ☐ 7 ☐ 8

Note: for each 13E400=* 2 domes and 13E401=* 3 domes 13Z16* are required.

Electrode domes for 13E400/401

Art. no.	Height	Qty.
☐ 13Z161	low	
☐ 13Z162	medium	
☐ 13Z163	high	

ErgoArm Hybrid Plus

Art. no.	Diameter	Color
☐ 12K44=45	45 mm	☐ (4) skin ☐ (7) black
☐ 12K44=50	50 mm	☐ (4) skin ☐ (7) black ☐ (11) brown* ☐ (15) dark brown*

ErgoArm Electronic Plus

Art. no.	Diameter	Color
☐ 12K50=45	45 mm	☐ (4) skin ☐ (7) black
☐ 12K50=50	50 mm	☐ (4) skin ☐ (7) black ☐ (11) brown* ☐ (15) dark brown*

DynamicArm

Art. no.	Diameter	Color
☐ 12K100N=45	45 mm	☐ (4) skin ☐ (7) black
☐ 12K100N=50	50 mm	☐ (4) skin ☐ (7) black ☐ (11) brown* ☐ (15) dark brown*

* Only for elbows with diameter 50 mm.

Rotation options with EQD connection

Description	Consists of
☐ Passive rotation	Coupling piece 10S4 and Coaxial plug 9E169
☐ Active rotation	Electric wrist rotator 10S17 and MyoRotronic 13E205
☐ Active rotation	Electric wrist rotator 10S17, when myosmart or DynamicArm is selected

Note: MyoRotronic is not shipped if 13E522 or 13E523=* is selected, except in combination with ErgoArm and "Active rotation".

Lamination ring

Art. no.	Diameter
☐ 10S1	☐ 45 mm ☐ 50 mm ☐ 54 mm

Energy-management and charger unit

Art. no. battery	Battery design	Art. no. w. charger	Color
☐ 757B35=4	MyoEnergy Integral	757L35	–
☐ 757B35=5	MyoEnergy Integral	757L35	–

Note: The charger will be selected automatically.

micelangelo and myosmart plus **for below & above elbow.** **Order form.**

Optional accessories

Skeo Up Daily Assist

Art. no.	Size	Qty.
☐ 14Y4	☐ 140 ☐ 160 ☐ 180 ☐ 200 ☐ 220	

Greifer

Art. no.	Qty.
☐ 8E33=9-1 DMC VariPlus	
☐ 8E34=9-1 DMC VariPlus	

Accessories socket

Art. no.	Description
☐ 12V10	Suction socket valve
☐ 99B13=21	PVC-pull-in tube
☐ 99B13=21-7	PVC-pull-in tube, black
☐ OC1560=TR	Easy Fit Arm for transradial fittings

Tool set for wrist assembly

Art. no.	Description
☐ 711M106=1	Mounting plate for wrist assembly
☐ 709H11=1	Spanner wrench for wrist assembly

Optional accessories elbow

Textile harness

Art. no.	Side	Size	Qty.
☐ 21A47=L-S	left	S	
☐ 21A47=L-M	left	M	
☐ 21A47=L-L	left	L	
☐ 21A47=R-S	right	S	
☐ 21A47=R-M	right	M	
☐ 21A47=R-L	right	L	
☐ OC1560=TH	Easy Fit Arm for transhumeral fittings.		

Accessoires myosmart

Art. no.	Description
☐ 13E101=1	Optional beeper for 13E522, 13E523=*
☐ 13E102=1	Elbow adapter cable for DynamicArm and ErgoArm 12K50

Optional accessories for gloves

Art. no.	Description
☐ 640F12	Special cleaner
☐ 640F13	Pump sprayer for 640F12 cleaner

USB-Adapter for charging (additional, for travelling)

☐ 757L43

Note: The 757L35 charger is not included in the scope of delivery of the 757L43 and needs to be ordered separately for a complete second charger with USB port must be ordered separately.

Myo evaluation tool

- ☐ 757M20-2=50 myosmart cuff
- ☐ 757M20-2=60 myosmart cuff

Accessoires DynamicArm

Art. no.	Description
☐ 646C42=V1.6	Elbowsoft 1.6
☐ 60X5	BionicLink
☐ 13E100	AnalogAdater for DynamicArm

Carbon Cover

Art. no.	Description
☐ 13Z157=45	for DynamicArm (diameter 45 mm)
☐ 13Z157=50	for DynamicArm (diameter 50 mm)

Comments:

Date _____ Place _____ Signature _____