

Liner portfolio

iconiq 3D printed silicone General information - Smooth and soft texture without stiffness - Textured inner surface for reliable suspension on the skin - Comfortable and skin-friendly - Easy to clean and maintain Indications - All residual limb types - Used with any socket	Skeo Silicone General information - Ideal for use in combination with a shuttle lock or lanyard system - Extremely durable - Easy to clean Indications - Residual limbs with a good soft tissue coverage - Use with a specific weight bearing socket (TT fitting)	Caleo Copolymer (TPE) General information - Best used in middle to higher activity with a cushion liner and a valve/sleeve combination, acceptable in lower activity with pin - Thermoformable, offers good protection against socket forces Indications - Many residual limb types; especially residual limbs with dry skin - Used with a total surface weight bearing socket (TT fitting)	Uneo Urethane General information - Best used in combination with valve or vacuum sockets - Material flow characteristics ensure even pressure distribution Indications - All residual limb types, especially sensitive, bony and/or scarred residual limbs - Used with a total surface weight bearing socket (TT fitting)
---	---	---	--

Amputation level	Liner family	Product image	Product name	Reference number	Characteristics	Distal connector	Textile cover	Sizes	Liner dimensions																																						
TT	iconiq		iconiq	6Y800 & 6Y801	- Uniform or tapered wall thickness - Available in three shore hardness levels for end-user optimized fit - Uniform compression - Textured inner surface for reliable suspension on the skin - Undercuts can be filled as required - In case of no outside textile cover: silky outer surface for easy donning and quick drying	—	● - OR —	Customized geometry based on the end-users' limb allow for individual fit and cover the full range of limb shapes, including challenging one's such as scar tissue indentations, undercuts, irregular contours etc.																																							
	Skeo		Skeo	6Y42	- Textile cover, easy to put on and take off without donning spray - Thin walls for flexibility - Reinforced to reduce distal elongation and positioning - Textured inner surface for reliable suspension on the skin	●	●	180 mm, 200 mm, 210 mm, 220 mm, 235 mm, 250 mm, 265 mm, 280 mm, 300 mm, 320 mm, 340 mm, 360 mm, 380 mm, 400 mm, 420 mm, 450 mm																																							
			Skeo Pure	6Y41	- Silky outer surface for easy donning and quick drying - Textured inner surface for reliable suspension on the skin - Distal matrix on locking version 6Y41 also available in 6 mm uniform thickness (e.g. 6Y41=280-6)	—	—	160 mm, 180 mm, 200 mm, 220 mm, 235 mm, 250 mm, 265 mm, 280 mm, 300 mm, 320 mm, 340 mm, 360 mm, 400 mm																																							
			Skeo	6Y43	- Reinforced to reduce distal elongation and positioning - Skinguard antibacterial (6Y75) - Textured inner surface for reliable suspension on the skin	●	—	120 mm, 140 mm, 160 mm, 180 mm, 200 mm, 210 mm, 220 mm, 235 mm, 250 mm, 265 mm, 280 mm, 300 mm, 320 mm, 340 mm, 360 mm, 380 mm, 400 mm, 420 mm, 450 mm																																							
				6Y70		●	●	160 mm, 180 mm, 200 mm, 210 mm, 220 mm, 235 mm, 250 mm, 265 mm, 280 mm, 300 mm, 320 mm, 340 mm, 360 mm, 380 mm, 400 mm																																							
			Skeo Skinguard	6Y75		●	●																																								
			Skeo Sealing TT	6Y112	- Integrated sealing ring for a long lifetime - Textured inner surface for reliable suspension on the skin - Silky outer surface for easy donning and quick drying - Orientation donning aid on the inner side	—	—	160x80 mm, 180x80 mm, 200x80 mm, 210x80 mm, 220x80 mm, 235x80 mm, 250x80 mm, 265x80 mm, 280x80 mm, 300x80 mm, 300x80 mm, 320x80 mm, 340x90 mm, 360x90 mm, 380x90 mm, 400x90 mm																																							
		Skeo Unique	6Y700	- Soft durometer silicone - Custom shaped from cast or scan for TT or TF - In case of no outside textile cover: silky outer surface for easy donning and quick drying - Easy to clean	● - OR —	● - OR —	Customized geometry based on limb. Ideal for challenging shapes and sizes, such as scarring, undercuts, knee flexion >15-35°, large circumference (>80cm) or long length (MPT to distal end >30cm) i.e. Syme, knee disarticulations, etc.	Customized																																							
	Caleo		Caleo	6Y90	- Thermoplastic elastomer (TPE) - Medical grade white oil 3mm thickness - Distal elongation reducing matrix (6Y90)	●	●	200 mm, 250 mm, 280 mm, 320 mm, 360 mm																																							
			Caleo	6Y92		—	●																																								
			Caleo 3D	6Y93=C6	- Thermoplastic elastomer (TPE) - Medical-grade white oil - Thinner wall thickness in rear proximal area (3mm) for easier knee flexion - Longer length 50 cm for ankle disarticulations (6Y93F=C6/L6*), also thicker cushion and locking (13 mm resp 21 mm)	—	●	<table><tr><th rowspan="2"></th><th colspan="2">4 cm from distal</th><th colspan="2">30 cm from distal</th></tr><tr><th>Min.</th><th>Max.</th><th>Min.</th><th>Max.</th></tr><tr><td>S</td><td>150</td><td>260</td><td>200</td><td>330</td></tr><tr><td>M</td><td>180</td><td>310</td><td>230</td><td>410</td></tr><tr><td>MP</td><td>200</td><td>310</td><td>330</td><td>500</td></tr><tr><td>L</td><td>230</td><td>350</td><td>360</td><td>560</td></tr><tr><td>LP</td><td>280</td><td>430</td><td>400</td><td>630</td></tr><tr><td>XL</td><td>330</td><td>500</td><td>430</td><td>690</td></tr></table>		4 cm from distal		30 cm from distal		Min.	Max.	Min.	Max.	S	150	260	200	330	M	180	310	230	410	MP	200	310	330	500	L	230	350	360	560	LP	280	430	400	630	XL	330	500	430	690
		4 cm from distal		30 cm from distal																																											
		Min.		Max.		Min.	Max.																																								
S		150	260	200	330																																										
M	180	310	230	410																																											
MP	200	310	330	500																																											
L	230	350	360	560																																											
LP	280	430	400	630																																											
XL	330	500	430	690																																											
	6Y93=L6	●	●																																												
	6Y93F=L6	●	●																																												
	Caleo 3D	6Y95=C6	- Thermoplastic elastomer (TPE) - Medical grade white oil - Thicker wall (6 mm) in front for better protection, and thinner wall (3 mm) rear - Flex-Zone in anterior material cover for optimal knee flexion ability - Improved gel formula combines durability and comfort, and can also be thermoformed - Soft more flexible 3D umbrella for locking liners and no umbrella for cushion liners	—	●	<table><tr><th rowspan="2"></th><th colspan="2">4 cm from distal</th><th colspan="2">30 cm from distal</th></tr><tr><th>Min.</th><th>Max.</th><th>Min.</th><th>Max.</th></tr><tr><td>S</td><td>150</td><td>260</td><td>200</td><td>330</td></tr><tr><td>M</td><td>180</td><td>310</td><td>230</td><td>410</td></tr><tr><td>MP</td><td>200</td><td>310</td><td>330</td><td>500</td></tr><tr><td>L</td><td>230</td><td>350</td><td>360</td><td>560</td></tr><tr><td>LP</td><td>280</td><td>430</td><td>400</td><td>630</td></tr><tr><td>XL</td><td>330</td><td>500</td><td>430</td><td>690</td></tr></table>		4 cm from distal		30 cm from distal		Min.	Max.	Min.	Max.	S	150	260	200	330	M	180	310	230	410	MP	200	310	330	500	L	230	350	360	560	LP	280	430	400	630	XL	330	500	430	690		
	4 cm from distal		30 cm from distal																																												
	Min.	Max.	Min.	Max.																																											
S	150	260	200	330																																											
M	180	310	230	410																																											
MP	200	310	330	500																																											
L	230	350	360	560																																											
LP	280	430	400	630																																											
XL	330	500	430	690																																											
	6Y95=L6	●	●																																												
	DVS Liner	6Y94a=*-F	- DVS Magnet connector - Soft TPE (Thermoplastic Elastomer) gel layer on the inside with partially exposed silicone layer on the outside for improved sealing and durability - Thinner wall thickness in rear proximal area (3mm) for optimal knee flexion.	●	●	<table><tr><th rowspan="2"></th><th colspan="2">4 cm from distal</th><th colspan="2">30 cm from distal</th></tr><tr><th>Min.</th><th>Max.</th><th>Min.</th><th>Max.</th></tr><tr><td>S</td><td>150</td><td>260</td><td>200</td><td>330</td></tr><tr><td>M</td><td>180</td><td>310</td><td>230</td><td>410</td></tr><tr><td>MP</td><td>200</td><td>310</td><td>330</td><td>500</td></tr><tr><td>L</td><td>230</td><td>350</td><td>360</td><td>560</td></tr><tr><td>LP</td><td>280</td><td>430</td><td>400</td><td>630</td></tr><tr><td>XL</td><td>330</td><td>500</td><td>430</td><td>690</td></tr></table>		4 cm from distal		30 cm from distal		Min.	Max.	Min.	Max.	S	150	260	200	330	M	180	310	230	410	MP	200	310	330	500	L	230	350	360	560	LP	280	430	400	630	XL	330	500	430	690		
	4 cm from distal		30 cm from distal																																												
	Min.	Max.	Min.	Max.																																											
S	150	260	200	330																																											
M	180	310	230	410																																											
MP	200	310	330	500																																											
L	230	350	360	560																																											
LP	280	430	400	630																																											
XL	330	500	430	690																																											
Uneo		Uneo	6Y520	- Ideal pressure distribution - Uniform wall thickness - In case of no outside textile cover: silky outer surface for easy donning and quick drying	—	● - OR —	190 mm, 230 mm, 280 mm, 330 mm, 400 mm																																								
		Uneo	6Y522	- Ideal pressure distribution - Tapered wall thickness - Skinguard Antibacterial option available (G) - Silky outer surface for easy donning and quick drying (only 6Y522)	—	—	190 mm, 210 mm, 230 mm, 250 mm, 280 mm, 310 mm																																								
			6Y523		—	●																																									
		Uneo 3D	6Y512	- Ideal pressure distribution - PreFixed 20 Degrees for easy bending - Anatomic specific geometry thicker in bony areas and thinner where flexion is needed - Skinguard Antibacterial option available (G) - Fresh version option available (P) - In case of no outside textile cover: silky outer surface for easy donning and quick drying	—	● - OR —	<table><tr><th rowspan="2">Circumference 4 cm from distal (g)</th><th colspan="3">Residual limb length from knee center (g)</th></tr><tr><th>75 mm</th><th>125 mm</th><th>175 mm</th></tr><tr><td>210 mm</td><td></td><td></td><td></td></tr><tr><td>235 mm</td><td></td><td></td><td></td></tr><tr><td>260 mm</td><td>*</td><td>*</td><td>*</td></tr><tr><td>285 mm</td><td>*</td><td>*</td><td>*</td></tr><tr><td>310 mm</td><td>*</td><td>*</td><td>*</td></tr><tr><td>335 mm</td><td>*</td><td>*</td><td>*</td></tr></table>	Circumference 4 cm from distal (g)	Residual limb length from knee center (g)			75 mm	125 mm	175 mm	210 mm				235 mm				260 mm	*	*	*	285 mm	*	*	*	310 mm	*	*	*	335 mm	*	*	*									
	Circumference 4 cm from distal (g)	Residual limb length from knee center (g)																																													
75 mm		125 mm	175 mm																																												
210 mm																																															
235 mm																																															
260 mm	*	*	*																																												
285 mm	*	*	*																																												
310 mm	*	*	*																																												
335 mm	*	*	*																																												
	Uneo Unique	6Y400	- Ideal pressure distribution - Custom shaped from cast or scan for TT or TF - Option for Fresh scent (P) or Skinguard Antibacterial (G) - In case of no outside textile cover: silky outer surface for easy donning and quick drying	● - OR —	● - OR —	Customized geometry based on limb. Ideal for challenging shapes and sizes, such as scarring, undercuts, knee flexion >15-35°, large circumference (>80cm) or long length (MPT to distal end >30cm) i.e. Syme, knee disarticulations, etc.	Customized																																								
TF	iconiq		iconiq	6Y830 & 6Y831	- Uniform or tapered wall thickness - Available in three shore hardness levels - Uniform compression - Textured inner surface for reliable suspension on the skin - Undercuts can be filled as - In case of no outside textile cover: silky outer surface for easy donning and quick drying	—	● - OR —	Customized geometry based on the end-users' limb allow for individual fit and cover the full range of limb shapes, including challenging one's such as scar tissue indentations, undercuts, irregular contours etc.																																							
	Skeo		Skeo	6Y80	- High circumferential elasticity - Integrated matrix for low lengthwise stretch - Skinguard Antibacterial (6Y85) - Textured inner surface for reliable suspension on the skin - Strain indicator on the outside of the liner for clear decision about size and shape of the liner	●	●	280 mm, 300 mm, 320 mm, 340 mm, 360 mm, 380 mm, 400 mm, 420 mm, 450 mm, 500 mm, 550 mm																																							
			Skeo Skinguard	6Y85		●	●																																								
			Skeo 3D	6Y87	3D shape, conical distal and standard proximal - Integrated Matrix for low lengthwise stretch - Textured inner surface for reliable suspension on the skin - Skinguard Antibacterial (6Y88) - Strain indicator on the outside of the liner for clear decision about size and shape of the liner	●	●	250 mm, 265 mm, 280 mm, 300 mm, 320 mm, 340 mm, 360 mm, 380 mm, 400 mm, 450 mm																																							
			Skeo 3D Skinguard	6Y88		●	●																																								
			Skeo Sealing	6Y110	- Integrated sealing ring for a long lifetime (10 cm or 17 cm) - Textured inner surface for reliable suspension on the skin - Silky outer surface for easy donning and quick drying - Orientation donning aid on the inner side 3D version is conical - Strain indicator on the outside of the liner for clear decision about size and shape of the liner - Selection chart is provided to simplify the decision about liner shape	—	—	<table><tr><th colspan="2">Ring height 10 cm</th><th colspan="2">Ring height 17 cm</th></tr><tr><th></th><th></th><th></th><th></th></tr><tr><td>250 x 10 mm (mm)</td><td>360 x 10</td><td>250 x 17 mm (mm)</td><td>360 x 17</td></tr><tr><td>265 x 10 mm (mm)</td><td>380 x 10</td><td>265 x 17 mm (mm)</td><td>380 x 17</td></tr><tr><td>280 x 10</td><td>400 x 10</td><td>280 x 17</td><td>400 x 17</td></tr><tr><td>300 x 10</td><td>450 x 10</td><td>300 x 17</td><td>450 x 17</td></tr><tr><td>320 x 10</td><td>500 x 10 mm (mm)</td><td>320 x 17</td><td>500 x 17 mm (mm)</td></tr><tr><td>340 x 10</td><td>550 x 10 mm (mm)</td><td>340 x 17</td><td>550 x 17 mm (mm)</td></tr></table>	Ring height 10 cm		Ring height 17 cm						250 x 10 mm (mm)	360 x 10	250 x 17 mm (mm)	360 x 17	265 x 10 mm (mm)	380 x 10	265 x 17 mm (mm)	380 x 17	280 x 10	400 x 10	280 x 17	400 x 17	300 x 10	450 x 10	300 x 17	450 x 17	320 x 10	500 x 10 mm (mm)	320 x 17	500 x 17 mm (mm)	340 x 10	550 x 10 mm (mm)	340 x 17	550 x 17 mm (mm)							
		Ring height 10 cm		Ring height 17 cm																																											
250 x 10 mm (mm)	360 x 10	250 x 17 mm (mm)	360 x 17																																												
265 x 10 mm (mm)	380 x 10	265 x 17 mm (mm)	380 x 17																																												
280 x 10	400 x 10	280 x 17	400 x 17																																												
300 x 10	450 x 10	300 x 17	450 x 17																																												
320 x 10	500 x 10 mm (mm)	320 x 17	500 x 17 mm (mm)																																												
340 x 10	550 x 10 mm (mm)	340 x 17	550 x 17 mm (mm)																																												
	Skeo Sealing 3D	6Y111	—	—																																											
	ProSeal Liner	6Y81 6Y81a=*-10	- High circumferential elasticity - Silky outer surface for easy donning and quick drying - Textured inner surface for reliable suspension on the skin - Also available with a 10 mm distal cushion (6Y81a=*-10)	—	—	280 mm, 300 mm, 320 mm, 340 mm, 360 mm, 380 mm, 400 mm, 420 mm, 450 mm, 500 mm, 550 mm																																									
	Skeo Unique	6Y730	- Soft durometer silicone - Custom shaped from cast or scan for TT or TF - In case of no outside textile cover: silky outer surface for easy donning and quick drying - Easy to clean	● - OR —	● - OR —	Customized geometry based on limb. Ideal for challenging shapes and sizes, such as scarring, undercuts, knee flexion >15-35°, large circumference (>80cm) or long length (MPT to distal end >30cm) i.e. Syme, knee disarticulations, etc.	Customized																																								
	Uneo Unique	6Y430	- Ideal pressure distribution - Custom shaped from cast or scan for TT or TF - Option for Fresh scent (P) or Skinguard Antibacterial (G) - In case of no outside textile cover: silky outer surface for easy donning and quick drying	● - OR —	● - OR —	Customized geometry based on limb. Ideal for challenging shapes and sizes, such as scarring, undercuts, knee flexion >15-35°, large circumference (>80cm) or long length (MPT to distal end >30cm) i.e. Syme, knee disarticulations, etc.	Customized																																								
TT TF TH TR	Compression liner		Compression Liner	6Y44	- Liner for compression therapy - Silky outer surface for easy donning and quick drying - Textured inner surface for reliable suspension on the skin - For use at different patients when being sterilised	—	—	160 mm, 180 mm, 200 mm, 220 mm, 235 mm, 250 mm, 265 mm, 280 mm, 300 mm, 320 mm, 340 mm, 360 mm, 380 mm, 400 mm, 450 mm, 500 mm, 550 mm																																							
			Compression Liner	6Y45	- Liner for compression therapy - Silky outer surface for easy donning and quick drying - Textured inner surface for reliable suspension on the skin - Liner can be used for the first fitting - For use at different patients when being sterilised	●	—	120 mm, 140 mm, 160 mm, 180 mm, 200 mm, 210 mm, 220 mm, 235 mm, 250 mm, 265 mm, 280 mm, 300 mm, 320 mm, 340 mm, 360 mm, 380 mm, 400 mm, 420 mm, 450 mm																																							