Doctor scan patient

Optional: Margin Marketing with iTero solutions, doctors can "make the margin" on the scan they send to you - resulting in more efficiency as you begin your workflow.

Optional

Return scans to doctor. If you find issues, like unclear margins or missing scans, you can immediately send scans back to the doctor with notes. An alert will appear on the iTero scanner so they can address the issue while the patient is still in the chair.

The iTero-enabled Workflow

MyiTero's single-uer interface allows you to direct your workflow in multiple ways, based on scan quality, case type, or production process.

Scan a transmitted in minutes

Lab receives scans

Once you recieve a scan, you control which step to take next. Need a model producer? Want to export the file directly to your preferred CAD software? You're in charge.

Optional iTero Services

Scan enhancement: There may be cases where a rescan isn't necessary but the scan could benefit from some enhancing (e.g. a margin could be more clearly defined) Just send the scan to Align for scan enhancement.

Precision milled models: If you're running at capacity and don't have time to 3D-print your own models, you can order a precise, high-quality milled model from Align.

