







flow, stop and go®



# NEOPERL® NEOSTRAHL CLINIC AERATORS

## STD size, laminar stream

### Models with restrictor

Color	Flow rate	STD insert	Female M22×1	Male M24×1
	Z 7.5 -9.0 l/min at 3 bar	<u>024063000</u>	<u>014098000</u>	<u>014198000</u>
	Z 7.5 -9.0 l/min at 3 bar	<u>024062000</u>	<u>014097000</u>	<u>014197000</u>
	Z 7.5 -9.0 l/min at 3 bar	<u>024060000</u>	<u>014095000</u>	<u>014195000</u>
	Z 7.5 -9.0 l/min at 3 bar	<u>024061000</u>	<u>014096000</u>	<u>014196000</u>