

ADVANCED WOUND CARE AT YOUR FINGERTIPS



RECREATING THE NATURAL WOUND HEALING ENVIRONMENT

ACTIGRAFT: THERE IS NO OTHER SUBSTITUTE

INTRODUCING ACTIGRAFT

The world's first innovative wound care product that enables health care providers to produce real-time whole blood clots to heal chronic ulcers.

COST-EFFECTIVE:

The only product that covers the wound for long periods of time without any care between applications.

EFFICIENT:

Provides the wound with an optimal, natural healing environment allowing the healing process to proceed unhindered.

SIMPLE:

Designed with nurses in mind, ActiGraft can be easily prepared at the patient's bedside in 12 minutes and requires no capital equipment.

AUTOGRAFT:

Created from the patient's blood, assuring safe and effective management of the wound.

PATIENT COMFORT:

Minimal handling of the wound.

ActiGraft is suitable for management of exuding cutaneous wounds, such as leg ulcers, pressure ulcers, diabetic ulcers, and mechanically or surgically-debrided wounds.



Scan here to view the ActiGraft procedure video

The American Medical Association (AMA) CPT Network affirmed in March 2019 that ActiGraft is described by existing skin substitute CPT codes, and is therefore eligible for coding and payment from CMS and other payers as a skin substitute.

RedDress offers innovative cost-effective wound-care solutions.

ActiGraft delivers the optimal solution for management of wounds in long term care facilities.



ACTIGRAFT PRODUCT



MICROSCOPIC VIEW OF ACTIGRAFT



reddressmedical.com

ACTIGRAFT: THERE IS NO SUBSTITUTE

INTRODUCING ACTIGRAFT

The world's first innovative wound care product that enables health care providers to produce real-time whole blood clots to heal chronic ulcers.

COST-EFFECTIVE:

The only product that covers the wound for long periods of time without any care between applications.

EFFICIENT:

Provides the wound with an optimal, natural healing environment allowing the healing process to proceed unhindered.

SIMPLE:

Designed with nurses in mind, ActiGraft can be easily prepared at the patient's bedside in 12 minutes and requires no capital equipment.

AUTOGRAFT:

Created from the patient's blood, assuring safe and effective management of the wound.

PATIENT COMFORT:

Minimal handling of the wound.

ActiGraft is suitable for management of exuding cutaneous wounds, such as leg ulcers, pressure ulcers, diabetic ulcers, and mechanically or surgically-debrided wounds.



Scan here to view the ActiGraft procedure video

The American Medical Association (AMA) CPT Network affirmed in March 2019 that ActiGraft is described by existing skin substitute CPT codes, and is therefore eligible for coding and payment from CMS and other payers as a skin substitute.

RedDress offers innovative cost-effective wound-care solutions. ActiGraft delivers the optimal solution for management of wounds in long term care facilities.

RedDress USA 822 A1A North, Ste 310 Ponte Vedra Beach, FL 32082 Tel.: +1.800.674.9615 Fax: +1.800.561.8127



ACTIGRAFT PRODUCT TOTAL CM² OF ACTIGRAFT: 28 CM²



MICROSCOPIC VIEW OF ACTIGRAFT

