

Embryon[®] BulbTip[™] Embryo Transfer Cath



For the introduction of human embryos into the uterine cavity during assisted reproduction

Ribbed Finger Guide providing grip and control ensuring accurate placing of the catheter tip in the uterine cavity

Depth insertion bead indicator in handle

- Bulb tip outer sheath with smooth bulb formed tip to guide the ultrafine inner catheter into the uterine cavity
- Ultrafine inner transfer catheter to minimise trauma and media volume on transfer
- Externally supported to ease loading and to aid insertion into the sheath during transfer.
- Sliding silicone marker to allow insertion depth to be pre-set prior to transfer
- Trial transfer catheters available
- Available in 18cm and 23cm
- Toxicity testing – MEA/LAL tested



DESCRIPTION	ORDER CODE	Length (cm)	Quantity	Outer Catheter OD (mm)	Inner Catheter		Volume mm ³
					OD (mm)	ID (mm)	
Embryon [®] BulbTip [™] Embryo Transfer Catheter	R57635-00-18	18	10	2.3	1.0	0.50	49.87
Embryon [®] BulbTip [™] Embryo Transfer Catheter	R57635-00-23	23	10	2.3	1.0	0.50	54.80
Embryon [®] BulbTip [™] Trial Transfer Catheter	R57636-00-18	18	10	2.3	1.0	0.5	49.87
Embryon [®] BulbTip [™] Trial Transfer Catheter	R57636-00-23	23	10	2.3	1.0	0.5	53.8