Do you know how to give your pet

CPR?



Confirm there's no heartbeat

Extend your pet's neck and gently clear any obstruction. Feel for the heartbeat by placing the hands on both sides of the chest behind the armpit/bent elbow for about 10 seconds.

Begin chest compressions

If your pet isn't breathing and no heartbeat is detected, start compressions. Place the pet on his/her side on a hard surface (ideally one that can also be used as a stretcher). Interlock hands with palms facing downward and compress the chest about half way for 60 compressions at a rate of 2 compressions per second.

Note: • Place hands over the heart (behind armpit/bent elbow) for pets less than 30 pounds.

 Place hands at the widest portion of the chest for pets greater than 30 pounds.

Start mouth-to-muzzle breathing

After 60 compressions, wrap your hands around the pet's muzzle to prevent air from escaping. Give 2 quick breaths, allowing the chest to expand and retract.

Resume compressions after the breaths and repeat this cycle until a heartbeat is detected or arrival to the veterinary facility. Check for a heartbeat after 4 cycles of CPR (approximately every 2 minutes).

Go immediately to your veterinarian clinic; call ahead so they're prepared for your arrival.

For more info on safety tips for pets, visit

MyPetHealthZone.com

